

# Modular Treatment Cabinet and Overhead Cabinet with Two Bottom Drawers and Doors

## Model 6117

The UMF Medical 6117 Treatment Cabinet has two lockable drawers and a two-door cabinet with a non-adjustable metal shelf in the lower section. The upper-section includes two lockable doors with glass panes, and one adjustable glass shelf. It comes equipped with a Kydex work surface and is the perfect complement to any exam room and treatment area.



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

### Key features at a glance

- › 20 gauge steel construction
- › One adjustable glass shelf in upper section
- › Removable, non-adjustable shelf in bottom section
- › Kydex top
- › Lock on upper doors (with 2 keys)
- › Heavy duty concealed door hinges
- › Safety release latch on upper left door
- › Adjustable leveling glides
- › 3-year limited warranty

### Specifications

|                            |                                 |
|----------------------------|---------------------------------|
| Width: 32" (81.3cm)        | Shipping Weight: 164lb (74.4kg) |
| Height: 65" (165.1cm)      | Shipping Dimensions:            |
| Depth: 16.25" (41.3cm)     | L35xW19xH72 (88.9cm x           |
| Net Weight: 134lb (60.8kg) | 48.3cm x 182.88cm)              |

### Accessories

502 - Double door lock (factory installed)  
 505 - Single drawer lock (factory installed)  
 507 - Additional keys (x2)