

Modular Treatment Cabinet with Two Bottom Drawers and Doors

Model 6120

UMF Medical's 6120 Modular Treatment Cabinet comes equipped with two bottom drawers, two-door cupboard with doors and non-adjustable shelf as well as a Kydex work surface. It is the perfect complement to any exam room and treatment area.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › 20 gauge steel construction
- › Kydex top with rail
- › Adjustable leveling glides
- › Removable, non-adjustable shelf
- › Heavy duty concealed door hinges
- › 3-year limited warranty

Specifications

Width: 27" (68.6cm)
Height: 33.5" (85.1cm)
Depth: 16.25" (41.3cm)
Net Weight: 73lb (33.1kg)

Shipping Weight: 87lb (39.5kg)
Shipping Dimensions:
L19xW30xH39 (48.3cm x
76.2cm x 90cm)

Accessories

502 - Double door lock (*factory installed*)
505 - Single drawer lock (*factory installed*)
507 - Additional keys (x2)