

Modular Treatment Cabinet with Four Drawers

Model 6204

UMF Medical's 6204 Treatment Cabinet has four drawers, a stainless steel push handle and twin wheel carpet casters. It comes equipped with a kydex work surface and is the perfect complement to any exam room and treatment area.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › 20 gauge steel construction
- › Stainless steel push handle
- › Kydex top with rail
- › 2" (5.1cm) twin wheel carpet casters
- › 3-year limited warranty

Specifications

Width: 25" (63.5cm)
Height: 34.25" (87cm)
Depth: 18" (45.7cm)
Net Weight: 89lb (40.4kg)

Shipping Weight: 145lb (65.8kg)
Shipping Dimensions:
L40xW48xH45 (101.6cm x
121.9cm x 114.3cm)

Accessories

122 – Locking swivel casters 2" (set of 4)
505 – Single drawer lock (factory installed)
507 – Additional keys (x2)