

Double Door, Double Lock Narcotic Cabinet with Two Shelves

Model 7781

UMF Medical's 7781 Narcotic Locker is manufactured with 20-gauge steel, offering double doors, double locks and two shelves.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › 20 gauge steel cabinet construction
- › 20 gauge double wall outer door
- › 14 gauge single wall inner door
- › Cam-style locks with 2 keys per lock
- › Heavy duty piano style door hinge
- › Two adjustable metal shelves
- › 3-year limited warranty

Specifications

Width: 12.125" (30.8cm)
Height: 15" (38.1cm)
Depth: 8.125" (20.6cm)
Net Weight: 17lb (7.7kg)

Shipping Weight: 22lb (9.8kg)
Shipping Dimensions:
L18xW15xH12 (41.7cm x
38.1cm x 30.5cm)

Accessories

506 – Master key feature (locks keyed differently, with master key)
507 – Additional keys (x2)
530 – 110V warning light (factory installed)
540 – 220V warning light (factory installed)