

Single Door, Double Lock Narcotic Cabinet with One Shelf

Model 7787

UMF Medical's 7787 Narcotic Locker is manufactured with 20-gauge steel, offering a single door, double lock and one shelf.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › 20 gauge steel cabinet construction
- › 20 gauge double wall outer door
- › Cam-style locks with 2 keys per lock
- › One adjustable shelf
- › Heavy-duty piano style door hinge
- › 3-year limited warranty

Specifications

Width: 14.125" (35.9cm)
Height: 12" (30.5cm)
Depth: 8.125" (20.6cm)
Net Weight: 14lb (6.4kg)

Shipping Weight: 19lb (8.6kg)
Shipping Dimensions:
L18xW15xH12 (45.7cm x
38.1cm x 30.5cm)

Accessories

506 - Master key feature (locks keyed differently, with master key)
507 - Additional keys (x2)
530 - 110V warning light (factory installed)
540 - 220V warning light (factory installed)