

# 48-Quart Red Waste Receptacle

## Model 1475R

The UMF Medical 1475R Waste Receptacle is made of heavy-duty steel and features a foot pedal operated lid. It also includes a bag securing mechanism for quick and easy clean up as well as free bio-hazard and infectious waste adhesive labels for OSHA labeling requirements. Comes with a 1-year limited warranty.



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

### Key features at a glance

- › 48 Quart (45.4L) Capacity
- › Flame Retardant
- › Unitized all-steel construction
- › Baked epoxy paint finish for added durability against scratching and stains
- › Full width hinge
- › Heavy-duty foot pedal and linkage
- › Self-closing lid
- › Exposed edges hemmed for safety
- › Bag securing mechanism
- › Free biohazard and infectious waste adhesive labels for OSHA labeling requirements
- › 1-year limited warranty

### Specifications

Height: 23.5" (59.7cm)	Shipping Weight: 19lb (8.6kg)
Width: 13" (33cm)	Shipping Dimensions: L25xW17xH16
Depth: 14" (35.6cm)	(63.5cm x 43.2cm x 40.7cm)
Net Weight: 15lb (6.8kg)	Color: Red