


32-Quart Stainless Steel Waste Receptacle

Model SS1474

The UMF Medical SS1474 Stainless Steel Waste Receptacle is made of heavy-duty steel and features a foot pedal operated lid. It also includes a bag securing mechanism for quick and easy clean up as well as free bio-hazard and infectious waste adhesive labels for OSHA labeling requirements. Comes with a 1-year limited warranty.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › Flame retardant
- › Unitized stainless steel construction
- › Resists stains and corrosion
- › Heavy-duty foot pedal and linkage
- › Full width hinge
- › Self-closing lid
- › Exposed edges hemmed for safety
- › Bag securing mechanism
- › Free biohazard and infectious waste adhesive labels for OSHA labeling requirements
- › 1-year limited warranty

Specifications

Height: 21" (53.3cm)
Width: 11.75" (29.8cm)
Depth: 13" (33cm)