

24" Four-Drawer Base Cabinet

Model 6004W

UMF Medical's 6004W 24" Four-Drawer Base Cabinet with full-extension roller bearing drawer glides is the perfect complement to any exam room's storage space. Our Wooden Cabinetry specialists work closely with you to understand your room layout and dimensions in order to offer a stylish, functional look to your exam room. All UMF Medical Wooden Cabinetry comes standard with interior and exterior laminate and is moisture and heat resistant.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Colors

Countertop



Titanium EV



Gesso Tracery



Classic Linen



Bronze Legacy

Cabinetry



White Cypress



Fusion Maple



Wild Cherry



Asian Night

Key features at a glance

- › Interior and exterior laminate
- › Moisture and heat resistant
- › Full extension roller bearing drawer glides
- › 3-year limited warranty

Specifications

Base cabinet only (including drawers):

Width: 24" (61cm)
Height: 34.5" (87.6cm)
Depth: 18" (45.7cm)

Base cabinet with countertop:

Width: 24" (61cm)
Height: 35.875" (91.1cm) – does not include height of backsplash
Height: 39.75" (101cm) – includes height of backsplash
Depth: 19" (48.3cm)

Accessories

- 610 – 49" Countertop with Sink cutout, pairs with 6014, 6024
- 620 – 25" Countertop with Sink cutout
- 625 – 25" Countertop without Sink cutout
- 636 – Sink assembly with Bowl, Faucet, and Strainer