

FusionONE™

Model 3002-500-105
Model 3002

For examinations

Every FusionONE™ model features a 500lb patient weight capacity and is available with Classic or Premium upholstery. Experience enhanced comfort, elegance, and softness in 16 stylish colors. The Premium upholstery model will elevate the aesthetics and comfort of any space.



Model 3002-500-105
includes Premium upholstery
(shown with optional armrests)

Key features at a glance

- › Power hi-lo and manual backrest
- › One-function foot control with hi-lo
- › Single column electro-magnetic actuator
- › 18"-21" paper roll holder (53.3cm)
- › 500lb (226.8kg) patient weight capacity
- › Large, reversible side drawer and large front drawer
- › Drain pan: plastic 1.5' deep (3.81cm)
- › Removable, seamless upholstered top with PreFixx® Protective Finish
- › Heavy-gauge steel extension leg pull-out to accommodate most patient heights
- › Stirrups with four lateral positions, adjustable length
- › Durable, steel construction
- › Easy-clean powder coat enamel surfaces
- › 3-year limited warranty

FusionONE™

Model 3002-500-105
Model 3002

Specifications

- › Height: min. 19" (48.3 cm) at seat
max. 32" (81.3 cm)
- › Width: 28" (71.12 cm)
- › Length: 58"-76" with foot section
extended (147.3 cm – 193.0 cm)
- › Leg rest: 18" W x 18" D (45.7 cm W x
45.7 cm D)
- › Storage Drawer (Front): 18.0" W x
12.0" D x 3.0" H (46 cm x 30 cm x 8 cm)
- › Storage Drawer (Side): 17.5" W x
21.0" D x 5.0" H (44 cm x 53 cm x 13 cm)
- › Single column electro-magnetic actuator
- › International model: 3002-01
- › Electrical: 120-240 VAC, 50-60 HZ
- › Weight capacity: 500lb (226.8kg)
- › Net weight: 300lb (136kg)
- › Shipping weight: 350lb (159kg)
- › Shipping dimensions: 30" W x
64" L x 31" H (76 cm x 163 cm x 79 cm)

Accessories

- › 251 – Bierhoff Knee Crutches (pair)
- › 581 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket
capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket
capabilities for #48955 (patient right)



Model 3002 includes
Classic upholstery

Upholstery and Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2026 United Metal Fabricators Inc. dba UMF Medical 05.26 All Rights Reserved. Prefix® is a registered trademark of Omnova Solutions. FusionONE is a trademark of United Metal Fabricators. 5-1371LT Rev B.



Ph: +1 814.266.8726 | Fax: +1 814.266.1870
1316 Eisenhower Blvd., Johnstown, PA 15904 USA

UMFmedical.com