

# Fusion**ONE**<sup>™</sup>

**Model 3002-500-105**  
**Model 3002**

## For examinations

With the accessibility of a power table and the affordability of a box table, the FusionONE Exam Table offers the features healthcare providers need most.

This hi-lo exam table with stirrups enables simple ADA access and ease of examination for patients of all ages and abilities.



Model 3002-500-105  
includes Premium upholstery  
(shown with optional armrests)

## Key features at a glance

- › Power hi-lo and manual backrest
- › One-function foot control with hi-lo
- › Single column electro-magnetic actuator
- › 18"-21" paper roll holder (53.3cm)
- › 500lb (226.8kg) patient weight capacity
- › ADA Accessible 19" (48.3 cm)
- › Large, reversible side drawer and large front drawer
- › Drain pan: plastic 1.5' deep (3.81cm)
- › Removable, seamless upholstered top with PreFixx<sup>®</sup> Protective Finish
- › Heavy-gauge steel extension leg pull-out to accommodate most patient heights
- › Stirrups with four lateral positions, adjustable length
- › Durable, steel construction
- › Easy-clean powder coat enamel surfaces
- › 3-year limited warranty

# FusionONE™

**Model 3002-500-105**  
**Model 3002**

## Specifications

- › Height: min. 19" (48.3 cm) at seat  
max. 32" (81.3 cm)
- › Width: 28" (71.12 cm)
- › Length: 58"-76" with foot section  
extended (147.3 cm – 193.0 cm)
- › Leg rest: 18" W x 18" D (45.7 cm W x  
45.7 cm D)
- › Storage Drawer (Front): 18.0" W x  
12.0" D x 3.0" H (46 cm x 30 cm x 8 cm)
- › Storage Drawer (Side): 17.5" W x  
21.0" D x 5.0" H (44 cm x 53 cm x 13 cm)
- › Single column electro-magnetic actuator
- › International model: 3002-01
- › Electrical: 120-240 VAC, 50-60 HZ
- › Weight capacity: 500lb (226.8kg)
- › Net weight: 300lb (136kg)
- › Shipping weight: 350lb (159kg)
- › Shipping dimensions: 30" W x  
64" L x 31" H (76 cm x 163 cm x 79 cm)

## Accessories

- › 251 – Bierhoff Knee Crutches (pair)
- › 581 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket  
capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket  
capabilities for #48955 (patient right)



Model 3002 includes  
Classic upholstery

## Upholstery and Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2022 United Metal Fabricators Inc. dba UMF Medical 08.22 All Rights Reserved. Prefix® is a registered trademark of Omnova Solutions. FusionONE is a trademark of United Metal Fabricators. 5-1371LT Rev A.



Ph: +1 814.266.8726 | Fax: +1 814.266.1870  
1316 Eisenhower Blvd., Johnstown, PA 15904 USA

[UMFmedical.com](http://UMFmedical.com)