

## Fusion**ONE**<sup>TM</sup>

### Power Hi-Lo, Manual Back Exam Chair with Foot Control, Stirrups, Drain Pan, Drawer Warmer, Pelvic Tilt & Receptacle

#### For examinations

With the accessibility of a power table and the affordability of a box table, the FusionONE Exam Chair offers the features healthcare providers need most – all in one chair!

This fully equipped, hi-lo exam chair can accommodate multiple specialties. It enables simple ADA access and ease of examination for patients of all ages and abilities.



Model 3003

#### Key features at a glance

- › Power hi-lo and manual backrest
- › One-function foot control with hi-lo
- › Single column electro-magnetic actuator
- › 18"-21" paper roll holder (53.3cm)
- › 500lb (226.8kg) patient weight capacity
- › Low-access ADA height 17"-19" (43.2cm-48.3cm)
- › Pelvic tilt (7°)
- › Drain pan: plastic 1.5" deep (3.81cm)
- › Large, reversible side drawer and large front drawer
- › Removable, seamless upholstered top with PreFixx<sup>®</sup> Protective Finish
- › Durable, galvanized steel construction
- › Hospital grade electrical receptacle with child safety cover
- › Heavy-gauge steel extension leg pull-out to accommodate most patient heights
- › Stirrups with four lateral positions, adjustable length
- › Easy-clean powder coat enamel surfaces
- › 3-year limited warranty

## Specifications

- › Height: min. 19" (±.5") 48.3 cm (±1.27 cm) at seat max. 32" (81.3 cm)
- › Width: 28" (71.12 cm)
- › Length: 58"-76" with foot section extended (147.3 cm – 193.0 cm)
- › Leg rest: 18" W x 18" D (45.7 cm W x 45.7 cm D)
- › Storage Drawer (Front): 18.0" W x 12.0" D x 3.0" H (46 cm x 30 cm x 8 cm)
- › Storage Drawer (Side): 17.5" W x 21.0" D x 5.0" H (44 cm x 53 cm x 13 cm)
- › Single column electro-magnetic actuator
- › International model: 3003-01
- › Electrical: 120-240 VAC, 50-60 HZ
- › Weight capacity: 500lb (226.8kg)
- › Net weight: 300lb (136kg)
- › Shipping weight: 350lb (159kg)
- › Shipping dimensions: 30" W x 64" L x 31" H (76 cm x 163 cm x 79 cm)



Model 3003  
shown with  
options

## Accessories

- › 580 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket capabilities for #48955 (patient right)

## Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2020 United Metal Fabricators Inc. dba UMF Medical 06.20 All Rights Reserved. Prefix® is a registered trademark of Omnova Solutions. FusionONE is a trademark of United Metal Fabricators. 5-1372LT Rev A.