

## Thermal Imaging Scanner Stand

**Model 7501** — for 4.5"-7" devices

**Model 7502** — for 7"-11" devices

**Model 7503** — for 9.5"-14.5" devices

The universal Thermal Imaging Scanner Stand was designed to securely hold thermal imaging scanners for hands-free temperature screening, tablets and smartphones. Lightweight and adjustable height with sturdy base, the stand features 360-degree lockable tilt with two axes of 360-degree rotation for infinitely adjustable device positioning to allow accurate, reliable thermal data capture. The wide, easy-grip, rubber-padded jaws securely hold almost any device — even in a protective case. Available in three sizes.



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**



### Key features at a glance

- › 360-degree lockable tilt with two axes of 360-degree rotation for infinitely adjustable device positioning to allow accurate, reliable thermal data capture
- › Wire management clip for routing the charging cord
- › Soft, rubber-padded grips prevent damage to the device
- › Diagonally grips device to ensure access to ports
- › Sand-cast base with rubber padding for reduced floor vibrations is sturdy yet light enough for transport
- › Diecast zinc clutch for height adjustments from 33"-60"

### Specifications

- › Color: Black
- › Base: 10" (25.4cm) diameter
- › Adjustable Height: 36" (91.4cm) – 60" (152.4cm)
- › Weight: 7.2lbs (3.3kg)