



# Fusion**ONE**<sup>TM</sup>

Series

## Power Exam Tables

**Accessible. Affordable.**

ADA Access in Every Exam Room.



See our quick,  
one-minute  
video on  
FusionONE!

[CLICK HERE](#)




**1 in 4**

adults is living with a disability.<sup>[1]</sup>



**1 in 2**

over 65 have a disability.<sup>[2]</sup>



**3 out of 4**

are overweight or obese.<sup>[3]</sup>



The population over age 65 is projected to grow from 48 million today to

**98 million**

in 2060<sup>[4]</sup>.

# Healthcare is changing

## And traditional exam rooms are changing with it

Traditional box exam tables can be difficult to access, posing safety risks for patients and clinicians.

And budget constraints often make it difficult to place ADA-accessible power exam chairs in every exam room.

But not anymore!

With the ACCESSIBILITY of a power chair and the AFFORDABILITY of a box table, the FusionONE Power Exam table offers the features healthcare providers need most.

This hi-lo power exam table enables simple ADA access and ease of examination for patients of all ages and abilities.

The FusionONE Series of exam tables are available with a manual or power backrest and with our OneTouch Wheelbase™ System. With a variety of models and accessories, the FusionONE Series allows for easy healthcare facility standardization.



◀ Model 3503

[1] US CDC [2] Harris Family Center for Disability and Health Policy / Western University of the Health Sciences [3] Modern Healthcare [4] The future is not what it used to be: Thoughts on the shape of the next normal, Kevin Sneader and Shubham Singhal, McKinsey & Company April 2020DC.org



	3001 3001-500-105	3002 3002-500-105	3003 3003-500-105	3002-500-300 3002-500-305	3003-500-300 3003-500-305	3502 3502-500-105	3503 3503-500-105	3502-500-300 3502-500-305	3503-500-300 3503-500-305
<b>Efficiencies</b>									
Power hi-lo	✓	✓	✓	✓	✓	✓	✓	✓	✓
Power backrest						✓	✓	✓	✓
Manual backrest	✓	✓	✓	✓	✓				
Foot control	✓	✓	✓	✓	✓	✓	✓	✓	✓
Front storage drawer	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reversible side drawer	✓	✓	✓			✓	✓	✓	✓
Plastic drain pan	✓	✓	✓	✓	✓	✓	✓	✓	✓
Paper roll holder	✓	✓	✓	✓	✓	✓	✓	✓	✓
Foldaway stirrups		✓	✓	✓	✓	✓	✓	✓	✓
Pelvic tilt 7°			✓		✓		✓		✓
<b>Infection Prevention</b>									
OneTouch WheelBase® System				✓	✓			✓	✓
Seamless upholstery with PreFixx® Protective Finish (Premium upholstery available on models ending in 05)	✓	✓	✓	✓	✓	✓	✓	✓	✓
Easy-clean powder coated enamel surfaces	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>Patient Experience</b>									
Drawer warmer			✓		✓		✓		✓
<b>Accessibility</b>									
500lb (226.8kg) patient weight capacity	✓	✓	✓	✓	✓	✓	✓	✓	✓
ADA Accessible 19" (48.3 cm)	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>Electrical</b>									
Hospital grade electrical receptacle			✓		✓		✓		✓
UL60601-1 and CAN/CSA-22.2 No. 601.1 classification	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>Quality</b>									
3-year limited warranty	✓	✓	✓	✓	✓	✓	✓	✓	✓
Durable galvanized steel construction	✓	✓	✓	✓	✓	✓	✓	✓	✓

**Specifications**

- › Height: min. 19" (48.3 cm) at seat max. 32" (81.3 cm)
- › Width: 28" (71.12 cm)
- › Length: 58"-76" with foot section extended (147.3 cm – 193.0 cm)
- › Leg rest: 18" W x 18" D (45.7 cm W x 45.7 cm D)
- › Storage Drawer (Front): 18.0" W x 12.0" D x 3.0" H (46 cm x 30 cm x 8 cm)
- › Storage Drawer (Side): 17.5" W x 21.0" D x 5.0" H (44 cm x 53 cm x 13 cm) (model dependent)
- › Single column electro-magnetic actuator
- › Electrical: 120-240 VAC, 50-60 HZ for 3000 series & 120-220 VAC, 50-60Hz for 3500 series
- › Amp range: 5 to 10.5 amp, depending on model
- › Weight capacity: 500lb (226.8kg)
- › Net weight: 325 lbs. (147kg) for 3000 series & 350 lbs. (159kg) for 3500 series
- › Shipping weight: 365 lbs. (166kg) for 3000 series & 390 lbs. (177kg) for 3500 series
- › Shipping dimensions: 30" W x 64" L x 31" H (76.2 cm x 162.6 cm x 78.7 cm)

**Accessories**

- › 251 – Bierhoff Knee Crutches (pair)
- › 581 – Hand control for 3000 series
- › 582 – Hand control for 3500 series
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket capabilities for #48955 (patient right)



FusionONE™ ProGlide™ FusionONE+™ ProGlide™

See our quick, one-minute video on the FusionONE ProGlide Series!

[CLICK HERE](#)

# The OneTouch WheelBase® System from UMF Medical

The OneTouch WheelBase® System from UMF Medical makes moving our ProGlide™ Series of products as easy as pushing a pedal. The integrated central-locking wheelbase with an intuitive foot pedal enables staff to move the equipment easily and quickly for cleaning, then lock it solidly into place for exams.

Model 3503-500-305 shown with optional armrests



## Upholstery & Colors

Every FusionONE and FusionONE+ model is available with Classic or Premium upholstery. Experience enhanced comfort, elegance, and softness in 16 stylish colors. The Premium upholstery model will elevate the aesthetics and comfort of any space.

### Colors



Ph: +1 814.266.8726 | Fax: +1 814.266.1870  
1316 Eisenhower Blvd., Johnstown, PA 15904 USA

[UMFmedical.com](http://UMFmedical.com)