



Fusion**ONE**<sup>TM</sup>  
Series  
Power Exam Tables





**1 in 4**

adults is living with a disability.<sup>[1]</sup>



**1 in 2**

over 65 have a disability.<sup>[2]</sup>



**3 out of 4**

are overweight or obese.<sup>[3]</sup>



The population over age 65 is projected to grow from 48 million today to

**98 million**

in 2060<sup>4</sup>.

# Healthcare is changing

## And traditional exam rooms are changing with it

The FusionONE Power Exam Table offers the features healthcare providers need most.

The OneTouch WheelBase® System from UMF Medical makes moving our ProGlide™ Series of products as easy as pushing a pedal. The integrated central-locking wheelbase with an intuitive foot pedal enables staff to move the equipment easily and quickly for cleaning, then lock it solidly into place for exams.



◀ Model 3503



	3001 3001-500-105	3002 3002-500-105	3003 3003-500-105	3002-500-300 3002-500-305	3003-500-300 3003-500-305	3502 3502-500-105	3503 3503-500-105	3502-500-300 3502-500-305	3503-500-300 3503-500-305
<b>Efficiencies</b>									
Power hi-lo	✓	✓	✓	✓	✓	✓	✓	✓	✓
Power backrest						✓	✓	✓	✓
Manual backrest	✓	✓	✓	✓	✓				
Foot control	✓	✓	✓	✓	✓	✓	✓	✓	✓
Front storage drawer	✓	✓	✓	✓	✓	✓	✓	✓	✓
Reversible side drawer	✓	✓	✓			✓	✓		
Plastic drain pan	✓	✓	✓	✓	✓	✓	✓	✓	✓
Paper roll holder	✓	✓	✓	✓	✓	✓	✓	✓	✓
Foldaway stirrups		✓	✓	✓	✓	✓	✓	✓	✓
Pelvic tilt 7°			✓		✓		✓		✓
500lb (226.8kg) patient weight capacity	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>Infection Prevention</b>									
OneTouch WheelBase® System				✓	✓			✓	✓
Seamless upholstery with PreFixx® Protective Finish (Premium upholstery available on models ending in 05)	✓	✓	✓	✓	✓	✓	✓	✓	✓
Easy-clean powder coated enamel surfaces	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>Patient Experience</b>									
Drawer warmer			✓		✓		✓		✓
<b>Electrical</b>									
Hospital grade electrical receptacle			✓		✓		✓		✓
UL60601-1 and CAN/CSA-22.2 No. 601.1 classification	✓	✓	✓	✓	✓	✓	✓	✓	✓
<b>Quality</b>									
3-year limited warranty	✓	✓	✓	✓	✓	✓	✓	✓	✓
Durable galvanized steel construction	✓	✓	✓	✓	✓	✓	✓	✓	✓

## Specifications

- › Height: min. 19" (48.3 cm) at seat max. 32" (81.3 cm)
- › Width: 28" (71.12 cm)
- › Length: 58"-76" with foot section extended (147.3 cm – 193.0 cm)
- › Leg rest: 18" W x 18" D (45.7 cm W x 45.7 cm D)
- › Storage Drawer (Front): 18.0" W x 12.0" D x 3.0" H (46 cm x 30 cm x 8 cm)
- › Storage Drawer (Side): 17.5" W x 21.0" D x 5.0" H (44 cm x 53 cm x 13 cm) (model dependent)
- › Single column electro-magnetic actuator
- › Electrical: 120-240 VAC, 50-60 HZ for 3000 series & 120-220 VAC, 50-60Hz for 3500 series
- › Amp range: 5 to 10.5 amp, depending on model
- › Weight capacity: 500lb (226.8kg)
- › Net weight: 325 lbs. (147kg) for 3000 series & 350 lbs. (159kg) for 3500 series
- › Shipping weight: 365 lbs. (166kg) for 3000 series & 390 lbs. (177kg) for 3500 series
- › Shipping dimensions: 30" W x 64" L x 31" H (76.2 cm x 162.6 cm x 78.7 cm)

## Accessories

- › 251 – Bierhoff Knee Crutches (pair)
- › 581 – Hand control for 3000 series
- › 582 – Hand control for 3500 series
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket capabilities for #48955 (patient right)



Fusion**ONE**<sup>™</sup> ProGlide<sup>™</sup> Fusion**ONE+**<sup>™</sup> ProGlide<sup>™</sup>

# The OneTouch WheelBase<sup>®</sup> System from UMF Medical

The OneTouch WheelBase<sup>®</sup> System from UMF Medical makes moving our ProGlide<sup>™</sup> Series of products as easy as pushing a pedal. The integrated central-locking wheelbase with an intuitive foot pedal enables staff to move the equipment easily and quickly for cleaning, then lock it solidly into place for exams.



◀ Model 3503-500-305 shown with optional armrests

## Upholstery & Colors

Every FusionONE and FusionONE+ model is available with Classic or Premium upholstery. Experience enhanced comfort, elegance, and softness in 16 stylish colors. The Premium upholstery model will elevate the aesthetics and comfort of any space.

### Colors



Ph: +1 814.266.8726 | Fax: +1 814.266.1870  
1316 Eisenhower Blvd., Johnstown, PA 15904 USA

[UMFmedical.com](http://UMFmedical.com)

©2026 United Metal Fabricators Inc. dba UMF Medical 05.26 All Rights Reserved. 5-1378LT Rev B