

Fusion**ONE**TM



FusionONE

Power Exam Chair

Accessible. Affordable.

ADA Access in Every Exam Room.




See our quick,
one-minute
video on
FusionONE!

[CLICK HERE](#)






1 in 4
adults is living
with a disability.^[1]



1 in 2
over 65 have
a disability.^[2]



3 out of 4
are overweight
or obese.^[3]



The population over age
65 is projected to grow
from 48 million today to
98 million
in 2060^[4].

Healthcare is changing.

And traditional exam rooms are changing with it.

Traditional box exam tables can be difficult to access, posing safety risks for patients and clinicians.

And budget constraints often make it difficult to place ADA-accessible power exam chairs in every exam room.

But not anymore!

With the **ACCESSIBILITY** of a power chair and the **AFFORDABILITY** of a box table, the FusionONE Power Exam Chair offers the features healthcare providers need most – all in one table!

This hi-lo power exam chair enables simple ADA access and ease of examination for patients of all ages and abilities.



[1] US CDC [2] Harris Family Center for Disability and Health Policy / Western University of the Health Sciences [3] Modern Healthcare [4] The future is not what it used to be: Thoughts on the shape of the next normal, Kevin Sneider and Shubham Singhal, McKinsey & Company April 2020DC.org



	3001	3002	3003
Infection Prevention			
Removable, seamless upholstery with PreFixx® Protective Finish	✓	✓	✓
Easy-clean powder coated enamel surfaces	✓	✓	✓
Efficiencies			
Foot control	✓	✓	✓
Power hi-lo with manual back	✓	✓	✓
Pelvic tilt (7°)			✓
Plastic drain pan			✓
Large front storage drawer	✓	✓	✓
Large, reversible side storage drawer	✓	✓	✓
Foldaway stirrups with adjustable length		✓	✓
Patient Experience			
Drawer warmer			✓
Comfort front leg cushion	✓	✓	✓
Accessibility			
Low-access ADA height 17"-19" (43.2cm-48.3)	✓	✓	✓
500lb (226.8kg) patient weight capacity	✓	✓	✓
Electrical			
Hospital grade receptacle with child safety cover			✓
UL60601-1 and CAN/CSA-22.2 No. 601.1 classification	✓	✓	✓
Quality			
3-year limited warranty	✓	✓	✓
Durable galvanized steel construction	✓	✓	✓

Specifications

- › Height: min. 19" (±.5") 48.3 cm (±1.27 cm) at seat max. 32" (81.3 cm)
- › Width: 28" (71.12 cm)
- › Length: 58"-76" with foot section extended (147.3 cm – 193.0 cm)
- › Leg rest: 18" W x 18" D (45.7 cm W x 45.7 cm D)
- › Storage Drawer (Front): 18.0" W x 12.0" D x 3.0" H (46 cm x 30 cm x 8 cm)
- › Storage Drawer (Side): 17.5" W x 21.0" D x 5.0" H (44 cm x 53 cm x 13 cm)
- › Single column electro-magnetic actuator
- › International models: 3001-01, 3002-01 and 3003-01
- › Electrical: 120-240 VAC, 50-60 HZ
- › Weight capacity: 500lb (226.8kg)
- › Net weight: 300lb (136kg)
- › Shipping weight: 350lb (159kg)
- › Shipping dimensions: 30" W x 64" L x 31" H (76 cm x 163 cm x 79 cm)

Accessories

- › 581 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket capabilities for #48955 (patient right)



Call today to schedule a live, online demo!

See the FusionONE Power Exam Chair, discover the many practical features of the design, and ask how this exam chair can work for you.



See our quick, one-minute video on FusionONE!

[CLICK HERE](#)



FusionONE Power Exam Chair

Accessible. Affordable.

ADA Access in Every Exam Room.

Colors

