

# power4010

head-centric procedure chair

## Procedure Chair with OneTouch Patient Positioning® System

### For head-centric practices

The power4010 is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers, such as ENT, Dermatology, and Plastic Surgery professionals. Designed for head-centric procedures, the power4010 accommodates the needs of most patients.



Model 4010-650-100  
Shown with optional  
Self-Leveling Arms

**PROTECTS** > patients from infection > facilities from damage > staff from injury

### Key features at a glance

- > Chair designed for head-centric practices
- > Slim back cushion for accessibility
- > Double-articulating headrest with height adjustability and interchangeable headrest options
- > 650 lb. (294.84 kg) patient weight capacity
- > Low height 18.5" (±0.5") 46.9cm (±1.27cm)
- > 4 power functions: hi-lo, back, tilt, leg lift & extension
- > Hand control with OneTouch Patient Positioning® System, 2 programmable positions and auto-return
- > Power seat that tilts up
- > Power leg lift for patient comfort
- > Concealed 18" paper roll holder
- > Trendelenberg, Lithotomy & Semi-Fowler positions
- > Removable, seamless upholstery with Prefix® Protective Finish
- > Hospital-grade electrical receptacle with child safety cover (patient left)



# power4010

head-centric procedure chair



## Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm) at seat max. 37" at seat 93.9cm
- › Seat width: 28" (71.12cm)
- › Back width: 24.5" (62.23cm)
- › Length: 70" – 80" with foot section extended (177.8cm – 203.2cm)
- › Leg rest: 28" W x 13" D (71.12cm W x 33.02cm D)
- › Dual electro-mechanical actuator system for efficient speed
- › Electrical: 120 VAC, 60HZ, 5 AMP
- › International model: 4010-650-101: 220/240V, 50/60Hz
- › Net weight: 410lb (185kg)
- › Shipping weight: 450lb (204kg)
- › Shipping dimensions: L64 x W30 x H41 (163cm x 76cm x 104cm)



Shown with options (047 Self-Leveling Arms)

## Accessories

- › 047 – Adjustable / self-leveling arms (pair)
- › 050 – Removable U-shaped headrest
- › 051 – Removable rectangular headrest
- › 234 – 2 Hook IV Pole Kit (includes clamp)
- › 235 – 4 Hook IV Pole Kit (includes clamp)
- › 238 – Fixed armboard package (EA)
- › 239 – Articulating armboard package (EA)
- › 473 – Welch Allyn light bkt capabilities for #48955 (patient left)
- › 474 – Welch Allyn light bkt capabilities for #48955 (patient right)
- › 830 – Patient safety strap (factory installed only)
- › 910 – Easy power table installation kit

## Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2024 United Metal Fabricators Inc. dba UMF Medical 12.24 All Rights Reserved. OneTouch Patient Positioning® is a registered trademark of United Metal Fabricators Inc. 5-1325LT Rev A



Ph: +1 814.266.8726 | Fax: +1 814.266.1870  
1316 Eisenhower Blvd., Johnstown, PA 15904 USA

UMFmedical.com

©2024 United Metal Fabricators Inc. dba UMF Medical