

FusionONE+™

Model 3502-500-105
Model 3502

For examinations

Every FusionONE+™ model features a 500lb patient weight capacity and a power back; this hi-lo exam table can accommodate multiple specialties. Experience enhanced comfort, elegance, and softness in 16 stylish colors. The Premium upholstery model will elevate the aesthetics and comfort of any space.



Model 3502-500-105 includes Premium upholstery (shown with optional armrests)

PROTECTS > patients from infection > facilities from damage > staff from injury

Key features at a glance

- > Power hi-lo and power backrest
- > Two-function foot control
- > 18"-21" (45.7 cm-53.3 cm) paper roll holder
- > 500 lb. (226.8 kg) patient weight capacity
- > Large, reversible side drawer and large front drawer
- > Drain pan: plastic 1.5" deep (3.8 cm)
- > Removable, seamless upholstered top with PreFixx® Protective Finish
- > Durable, steel construction
- > Heavy-gauge steel pull-out leg rest
- > Foldaway stirrups with four lateral positions and adjustable length
- > Easy-clean powder coat enamel surfaces
- > 3-year limited warranty

FusionONE+™

Model 3502-500-105
Model 3502

Specifications

- › Height: min. 19" (48.3 cm) at seat max. 32" (81.3 cm)
- › Width: 28" (71.1 cm)
- › Length: 58"-76" (147.3 cm-193 cm) with leg rest extended
- › Leg rest: 18" W x 18" D (45.7 cm W x 45.7 cm D)
- › Storage Drawer (Front): 18" W x 12" D x 3" H (45.7 cm x 30.5 cm x 7.6 cm)
- › Storage Drawer (Side): 17.5" W x 21" D x 5" H (44.5 cm x 53.3 cm x 12.7 cm)
- › Single-column, electro-magnetic actuator
- › Electrical: 120 VAC 10 AMP / 220 VAC 5.5 AMP, 50-60 Hz
- › Weight capacity: 500 lbs. (226.8kg)
- › Net weight: 350 lbs. (159kg)
- › Shipping weight: 390 lbs. (177kg)
- › Shipping dimensions: 30" W x 64" L x 31" H (76.2 cm x 162.6 cm x 78.7 cm)

Accessories

- › 251 – Bierhoff Knee Crutches (pair)
- › 582 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket capabilities for #48955 (patient right)



Model 3502 includes Classic upholstery

Upholstery and Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2026 United Metal Fabricators Inc. dba UMF Medical 05.26 All Rights Reserved. PreFixx® is a registered trademark of Omnova Solutions. FusionONE is a trademark of United Metal Fabricators. 5-1372LT Rev B.



Ph: +1 814.266.8726 | Fax: +1 814.266.1870
1316 Eisenhower Blvd., Johnstown, PA 15904 USA

UMFmedical.com

©2026 United Metal Fabricators Inc. dba UMF Medical