

FusionONE+™

Model 3503-500-105
Model 3503

For examinations

Fusing the accessibility of a power table and the affordability of a fixed-height box table, the FusionONE+ Power Exam Table offers the features healthcare providers need most.

Featuring a power back, this hi-lo exam table can accommodate multiple specialties. It enables simple ADA access and ease of examination for patients of all ages and abilities. The power backrest makes it easy to lie patients back and sit them up safely and comfortably.



Model 3503-500-105 includes Premium upholstery (shown with optional armrests)

PROTECTS > patients from infection > facilities from damage > staff from injury

Key features at a glance

- > Power hi-lo and power backrest
- > Two-function foot control
- > 18"-21" (45.7 cm-53.3 cm) paper roll holder
- > 500 lb. (226.8 kg) patient weight capacity
- > ADA Accessible 19" (48.3 cm)
- > Pelvic tilt (7°)
- > Drain pan: plastic 1.5" deep (3.8 cm)
- > Large, reversible side drawer and large front drawer
- > Removable, seamless upholstered top with PreFixx® Protective Finish
- > Durable, steel construction
- > Electrical receptacle
- > Drawer warmer
- > Heavy-gauge steel pull-out leg extension
- > Foldaway stirrups with four lateral positions and adjustable length
- > Easy-clean powder coat enamel surfaces
- > 3-year limited warranty

FusionONE+™

Model 3503-500-105
Model 3503

Specifications

- › Height: min. 19" (48.3 cm) at seat max. 32" (81.3 cm)
- › Width: 28" (71.1 cm)
- › Length: 58"-76" (147.3 cm-193 cm) with leg rest extended
- › Leg rest: 18" W x 18" D (45.7 cm W x 45.7 cm D)
- › Storage Drawer (Front): 18" W x 12" D x 3" H (45.7 cm x 30.5 cm x 7.6 cm)
- › Storage Drawer (Side): 17.5" W x 21" D x 5" H (44.5 cm x 53.3 cm x 12.7 cm)
- › Single-column, electro-magnetic actuator
- › Electrical: 120 VAC 10.5 AMP / 220 VAC 5.8 AMP, 50-60Hz
- › Weight capacity: 500 lbs. (226.8kg)
- › Net weight: 350 lbs. (159kg)
- › Shipping weight: 390 lbs. (177kg)
- › Shipping dimensions: 30" W x 64" L x 31" H (76.2 cm x 162.6 cm x 78.7 cm)

Accessories

- › 251 – Bierhoff Knee Crutches (pair)
- › 582 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket capabilities for #48955 (patient right)



Model 3503 includes Classic upholstery

Upholstery and Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2023 United Metal Fabricators Inc. dba UMF Medical 07.23 All Rights Reserved. PreFixx® is a registered trademark of Omnova Solutions. FusionONE is a trademark of United Metal Fabricators. 5-1372LT Rev A.



Ph: +1 814.266.8726 | Fax: +1 814.266.1870
1316 Eisenhower Blvd., Johnstown, PA 15904 USA

UMFmedical.com

©2023 United Metal Fabricators Inc. dba UMF Medical