

# FusionONE+™

**Model 3503-500-105**  
**Model 3503**

## For examinations

Every FusionONE+™ model features a 500lb patient weight capacity and a power back; this hi-lo exam table can accommodate multiple specialties. Experience enhanced comfort, elegance, and softness in 16 stylish colors. The Premium upholstery model will elevate the aesthetics and comfort of any space.



Model 3503-500-105 includes Premium upholstery (shown with optional armrests)

**PROTECTS** > patients from infection > facilities from damage > staff from injury

## Key features at a glance

- > Power hi-lo and power backrest
- > Two-function foot control
- > 18"-21" (45.7 cm-53.3 cm) paper roll holder
- > 500 lb. (226.8 kg) patient weight capacity
- > Pelvic tilt (7°)
- > Drain pan: plastic 1.5" deep (3.8 cm)
- > Large, reversible side drawer and large front drawer
- > Removable, seamless upholstered top with PreFixx® Protective Finish
- > Durable, steel construction
- > Electrical receptacle
- > Drawer warmer
- > Heavy-gauge steel pull-out leg extension
- > Foldaway stirrups with four lateral positions and adjustable length
- > Easy-clean powder coat enamel surfaces
- > 3-year limited warranty

# FusionONE+™

**Model 3503-500-105**  
**Model 3503**

## Specifications

- › Height: min. 19" (48.3 cm) at seat max. 32" (81.3 cm)
- › Width: 28" (71.1 cm)
- › Length: 58"-76" (147.3 cm-193 cm) with leg rest extended
- › Leg rest: 18" W x 18" D (45.7 cm W x 45.7 cm D)
- › Storage Drawer (Front): 18" W x 12" D x 3" H (45.7 cm x 30.5 cm x 7.6 cm)
- › Storage Drawer (Side): 17.5" W x 21" D x 5" H (44.5 cm x 53.3 cm x 12.7 cm)
- › Single-column, electro-magnetic actuator
- › Electrical: 120 VAC 10.5 AMP / 220 VAC 5.8 AMP, 50-60Hz
- › Weight capacity: 500 lbs. (226.8kg)
- › Net weight: 350 lbs. (159kg)
- › Shipping weight: 390 lbs. (177kg)
- › Shipping dimensions: 30" W x 64" L x 31" H (76.2 cm x 162.6 cm x 78.7 cm)

## Accessories

- › 251 – Bierhoff Knee Crutches (pair)
- › 582 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket capabilities for #48955 (patient right)



Model 3503 includes Classic upholstery

## Upholstery and Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2026 United Metal Fabricators Inc. dba UMF Medical 05.26 All Rights Reserved. PreFixx® is a registered trademark of Omnova Solutions. FusionONE is a trademark of United Metal Fabricators. 5-1372LT Rev B.



Ph: +1 814.266.8726 | Fax: +1 814.266.1870  
1316 Eisenhower Blvd., Johnstown, PA 15904 USA

**UMFmedical.com**

©2026 United Metal Fabricators Inc. dba UMF Medical