

Stainless Steel Utility Table with One Drawer, One Shelf, and Three-sided Guard Rail

Model SS8153

UMF Medical Stainless Steel Utility Tables come with swivel casters, and have a 3-year limited warranty.

The SS8153 has a single drawer and storage shelf below the drawer. A 3-sided guard rail surrounds the top work surface.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › Welded stainless steel construction
- › Stainless steel drawer
- › Three-sided guard rail on top
- › Convenient storage shelf
- › 2" (5.1cm) swivel casters
- › 3-year limited warranty

Specifications

Width: 20" (50.8cm)
Height: 35.3" (89.5cm)
Depth: 16" (40.6cm)

Accessories

117 – Swivel casters locking 2" (set of 4)
165 – Stainless Steel push handle
480 – Rubber crutch tips in place of casters (table height will be 2" (5.1cm)) lower
505 – Single drawer lock (factory installed)