

FusionONE™

ProGlide™

Model 3003-500-305
Model 3003-500-300

For examinations

With the accessibility of a power table and the affordability of a box table, the FusionONE ProGlide offers the features healthcare providers need most.

This fully equipped, hi-lo exam table can accommodate multiple specialties. It enables simple ADA access and ease of examination for patients of all ages and abilities.

The OneTouch WheelBase® System from UMF Medical makes moving our ProGlide™ Series of products as easy as pushing a pedal. The integrated central-locking wheelbase with an intuitive foot pedal enables staff to move the equipment easily and quickly for cleaning, then lock it solidly into place for exams.



Model 3003-500-305 includes Premium upholstery (shown with optional armrests)

Key features at a glance

- › OneTouch WheelBase® System
- › Power hi-lo and manual backrest
- › One-function foot control with hi-lo
- › Single column electro-magnetic actuator
- › 18"-21" paper roll holder (53.3cm)
- › 500lb (226.8kg) patient weight capacity
- › ADA Accessible 19" (48.3 cm)
- › Pelvic tilt (7°)
- › Drain pan: plastic 1.5" deep (3.81cm)
- › Large front drawer
- › Removable, seamless upholstered top with PreFixx® Protective Finish
- › Durable, steel construction
- › Hospital grade electrical receptacle
- › Drawer warmer
- › Heavy-gauge steel extension leg pull-out to accommodate most patient heights
- › Stirrups with four lateral positions, adjustable length
- › Easy-clean powder coat enamel surfaces
- › 3-year limited warranty

FusionONE™

ProGlide™

Model 3003-500-305
Model 3003-500-300

Specifications

- › Height: min. 19" (48.3 cm)
at seat max. 32" (81.3 cm)
- › Width: 28" (71.12 cm)
- › Length: 58"-76" with foot section
extended (147.3 cm – 193.0 cm)
- › Leg rest: 18" W x 18" D (45.7 cm W
x 45.7 cm D)
- › Storage Drawer (Front): 18.0" W x
12.0" D x 3.0" H (46 cm x 30 cm x 8 cm)
- › Single column electro-magnetic actuator
- › International model: 3003-01
- › Electrical: 120-240 VAC, 50-60 HZ
- › Weight capacity: 500lb (226.8kg)
- › Net weight: 300lb (136kg)
- › Shipping weight: 350lb (159kg)
- › Shipping dimensions: 30" W x
64" L x 31" H (76 cm x 163 cm x 79 cm)

Accessories

- › 251 – Bierhoff Knee Crutches (pair)
- › 581 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket
capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket
capabilities for #48955 (patient right)



Model 3003-500-300
includes Classic upholstery

Upholstery and Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2023 United Metal Fabricators Inc. dba UMF Medical 07.23 All Rights Reserved. Prefix® is a registered trademark of Omnova Solutions. FusionONE is a trademark of United Metal Fabricators. 5-1372LT Rev A.



Ph: +1 814.266.8726 | Fax: +1 814.266.1870
 1316 Eisenhower Blvd., Johnstown, PA 15904 USA

UMFmedical.com