

# FusionONE™

---

## ProGlide™

**Model 3003-500-305**  
**Model 3003-500-300**

### For examinations

Every FusionONE™ model features a 500lb patient weight capacity and is available with Classic or Premium upholstery. Experience enhanced comfort, elegance, and softness in 16 stylish colors. The Premium upholstery model will elevate the aesthetics and comfort of any space.

The OneTouch WheelBase® System from UMF Medical makes moving our ProGlide™ Series of products as easy as pushing a pedal. The integrated central-locking wheelbase with an intuitive foot pedal enables staff to move the equipment easily and quickly for cleaning, then lock it solidly into place for exams.



Model 3003-500-305  
includes Premium upholstery  
(shown with optional armrests)

### Key features at a glance

- › OneTouch WheelBase® System
- › Power hi-lo and manual backrest
- › One-function foot control with hi-lo
- › Single column electro-magnetic actuator
- › 18"-21" paper roll holder (53.3cm)
- › 500lb (226.8kg) patient weight capacity
- › Pelvic tilt (7°)
- › Drain pan: plastic 1.5" deep (3.81cm)
- › Large front drawer
- › Removable, seamless upholstered top with PreFixx® Protective Finish
- › Durable, steel construction
- › Hospital grade electrical receptacle
- › Drawer warmer
- › Heavy-gauge steel extension leg pull-out to accommodate most patient heights
- › Stirrups with four lateral positions, adjustable length
- › Easy-clean powder coat enamel surfaces
- › 3-year limited warranty

# FusionONE™

ProGlide™

**Model 3003-500-305**  
**Model 3003-500-300**

## Specifications

- › Height: min. 19" (48.3 cm)  
at seat max. 32" (81.3 cm)
- › Width: 28" (71.12 cm)
- › Length: 58"-76" with foot section  
extended (147.3 cm – 193.0 cm)
- › Leg rest: 18" W x 18" D (45.7 cm W  
x 45.7 cm D)
- › Storage Drawer (Front): 18.0" W x  
12.0" D x 3.0" H (46 cm x 30 cm x 8 cm)
- › Single column electro-magnetic actuator
- › International model: 3003-01
- › Electrical: 120-240 VAC, 50-60 HZ
- › Weight capacity: 500lb (226.8kg)
- › Net weight: 300lb (136kg)
- › Shipping weight: 350lb (159kg)
- › Shipping dimensions: 30" W x  
64" L x 31" H (76 cm x 163 cm x 79 cm)

## Accessories

- › 251 – Bierhoff Knee Crutches (pair)
- › 581 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket  
capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket  
capabilities for #48955 (patient right)



Model 3003-500-300  
includes Classic upholstery

## Upholstery and Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2026 United Metal Fabricators Inc. dba UMF Medical 05.26 All Rights Reserved. Prefix® is a registered trademark of Omnova Solutions. FusionONE is a trademark of United Metal Fabricators. 5-1372LT Rev B.



Ph: +1 814.266.8726 | Fax: +1 814.266.1870  
 1316 Eisenhower Blvd., Johnstown, PA 15904 USA

**UMFmedical.com**