

FusionONE+™

ProGlide™

Model 3502-500-300

For examinations

Fusing the accessibility of a power table and the affordability of a fixed-height box table, the FusionONE+ ProGlide offers the features healthcare providers need most.

Featuring a power back, this hi-lo exam table can accommodate multiple specialties. It enables simple ADA access and ease of examination for patients of all ages and abilities. The power backrest makes it easy to lie patients back and sit them up safely and comfortably.

The OneTouch WheelBase® System from UMF Medical makes moving our ProGlide™ Series of products as easy as pushing a pedal. The integrated central-locking wheelbase with an intuitive foot pedal enables staff to move the equipment easily and quickly for cleaning, then lock it solidly into place for exams.



Key features at a glance

- › OneTouch WheelBase® System
- › ADA Accessible 19" (48.3 cm)
- › Heavy-gauge steel pull-out leg rest
- › Power hi-lo and power backrest
- › Large front drawer
- › Foldaway stirrups with four lateral positions and adjustable length
- › Two-function foot control
- › Drain pan: plastic 1.5" deep (3.8 cm)
- › Easy-clean powder coat enamel surfaces
- › 18"-21" (45.7 cm-53.3 cm) paper roll holder
- › Removable, seamless upholstered top with PreFixx® Protective Finish
- › 3-year limited warranty
- › 500 lb. (226.8 kg) patient weight capacity
- › Durable, galvanized steel construction

FusionONE+™

ProGlide™



Specifications

- › Height: min. 19" (48.3 cm) at seat max. 32" (81.3 cm)
- › Width: 28" (71.1 cm)
- › Length: 58"-76" (147.3 cm-193 cm) with leg rest extended
- › Leg rest: 18" W x 18" D (45.7 cm W x 45.7 cm D)
- › Storage Drawer (Front): 18" W x 12" D x 3" H (45.7 cm x 30.5 cm x 7.6 cm)
- › Single-column, electro-magnetic actuator
- › Electrical: 120 VAC 10 AMP / 220 VAC 5.5 AMP, 50-60 Hz
- › Weight capacity: 500 lbs. (226.8kg)
- › Net weight: 350 lbs. (159kg)
- › Shipping weight: 390 lbs. (177kg)
- › Shipping dimensions: 30" W x 64" L x 31" H (76.2 cm x 162.6 cm x 78.7 cm)



Model 3502-500-300 shown with options

Accessories

- › 251 – Bierhoff Knee Crutches (pair)
- › 582 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket capabilities for #48955 (patient right)

Upholstery and Colors

Every model is available with Classic or Premium upholstery and in 16 colors.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2023 United Metal Fabricators Inc. dba UMF Medical 07.23 All Rights Reserved. Prefix® is a registered trademark of Omnova Solutions. FusionONE is a trademark of United Metal Fabricators. 5-1372LT Rev A.



Ph: +1 814.266.8726 | Fax: +1 814.266.1870
1316 Eisenhower Blvd., Johnstown, PA 15904 USA

UMFmedical.com

©2023 United Metal Fabricators Inc. dba UMF Medical