

Fusion**ONE+**[™] ProGlide[™]

Model 3503-500-305

Model 3503-500-300

For examinations

Every FusionONE+[™] model features a 500lb patient weight capacity and is available with Classic or Premium upholstery. Experience enhanced comfort, elegance, and softness in 16 stylish colors. The Premium upholstery model will elevate the aesthetics and comfort of any space.

The OneTouch WheelBase[®] System from UMF Medical makes moving our ProGlide[™] Series of products as easy as pushing a pedal. The integrated central-locking wheelbase with an intuitive foot pedal enables staff to move the equipment easily and quickly for cleaning, then lock it solidly into place for exams.



Model 3503-500-305
includes Premium upholstery
(shown with optional armrests)

Key features at a glance

- › OneTouch WheelBase[®] System
- › Power hi-lo and power backrest
- › Two-function foot control
- › 18"-21" (45.7 cm-53.3 cm) paper roll holder
- › 500 lb. (226.8 kg) patient weight capacity
- › Pelvic tilt (7°)
- › Drain pan: plastic 1.5" deep (3.8 cm)
- › Large front drawer
- › Removable, seamless upholstered top with PreFixx[®] Protective Finish
- › Durable, steel construction
- › Drawer warmer
- › Electrical receptacle
- › Heavy-gauge steel pull-out leg extension
- › Foldaway stirrups with four lateral positions and adjustable length
- › Easy-clean powder coat enamel surfaces
- › 3-year limited warranty

FusionONE+™

ProGlide™

Model 3503-500-305
Model 3503-500-300

Specifications

- › Height: min. 19" (48.3 cm)
at seat max. 32" (81.3 cm)
- › Width: 28" (71.1 cm)
- › Length: 58"-76" (147.3 cm-193 cm)
with leg rest extended
- › Leg rest: 18" W x 18" D (45.7 cm W
x 45.7 cm D)
- › Storage Drawer (Front): 18" W x 12" D
x 3" H (45.7 cm x 30.5 cm x 7.6 cm)
- › Single-column, electro-magnetic actuator
- › Electrical: 120 VAC 10.5 AMP / 220 VAC
5.8 AMP, 50-60Hz
- › Weight capacity: 500 lbs. (226.8kg)
- › Net weight: 350 lbs. (159kg)
- › Shipping weight: 390 lbs. (177kg)
- › Shipping dimensions: 30" W x 64" L x 31" H
(76.2 cm x 162.6 cm x 78.7 cm)

Accessories

- › 251 – Bierhoff Knee Crutches (pair)
- › 582 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket
capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket
capabilities for #48955 (patient right)



Model 3503-500-300
includes Classic upholstery

Upholstery and Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2026 United Metal Fabricators Inc. dba UMF Medical 05.26 All Rights Reserved. PreFixx® is a registered trademark of Omnova Solutions. FusionONE is a trademark of United Metal Fabricators. 5-1372LT Rev B.



Ph: +1 814.266.8726 | Fax: +1 814.266.1870
 1316 Eisenhower Blvd., Johnstown, PA 15904 USA

UMFmedical.com

©2026 United Metal Fabricators Inc. dba UMF Medical