

# FusionONE+™ ProGlide™

**Model 3503-500-305**


**Model 3503-500-300**

## For examinations

Fusing the accessibility of a power table and the affordability of a fixed-height box table, the FusionONE+ ProGlide offers the features healthcare providers need most.

Featuring a power back, this hi-lo exam table can accommodate multiple specialties. It enables simple ADA access and ease of examination for patients of all ages and abilities. The power backrest makes it easy to lie patients back and sit them up safely and comfortably.

The OneTouch WheelBase® System from UMF Medical makes moving our ProGlide™ Series of products as easy as pushing a pedal. The integrated central-locking wheelbase with an intuitive foot pedal enables staff to move the equipment easily and quickly for cleaning, then lock it solidly into place for exams.



Model 3503-500-305 includes Premium upholstery (shown with optional armrests)

## Key features at a glance

- › OneTouch WheelBase® System
- › Power hi-lo and power backrest
- › Two-function foot control
- › 18"-21" (45.7 cm-53.3 cm) paper roll holder
- › 500 lb. (226.8 kg) patient weight capacity
- › ADA Accessible 19" (48.3 cm)
- › Pelvic tilt (7°)
- › Drain pan: plastic 1.5" deep (3.8 cm)
- › Large front drawer
- › Removable, seamless upholstered top with PreFixx® Protective Finish
- › Durable, steel construction
- › Electrical receptacle
- › Heavy-gauge steel pull-out leg extension
- › Foldaway stirrups with four lateral positions and adjustable length
- › Easy-clean powder coat enamel surfaces
- › 3-year limited warranty

# FusionONE+™

## ProGlide™

**Model 3503-500-305**

**Model 3503-500-300**

### Specifications

- › Height: min. 19" (48.3 cm)  
at seat max. 32" (81.3 cm)
- › Width: 28" (71.1 cm)
- › Length: 58"-76" (147.3 cm-193 cm)  
with leg rest extended
- › Leg rest: 18" W x 18" D (45.7 cm W  
x 45.7 cm D)
- › Storage Drawer (Front): 18" W x 12" D  
x 3" H (45.7 cm x 30.5 cm x 7.6 cm)
- › Single-column, electro-magnetic actuator
- › Electrical: 120 VAC 10.5 AMP / 220 VAC  
5.8 AMP, 50-60Hz
- › Weight capacity: 500 lbs. (226.8kg)
- › Net weight: 350 lbs. (159kg)
- › Shipping weight: 390 lbs. (177kg)
- › Shipping dimensions: 30" W x 64" L x 31" H  
(76.2 cm x 162.6 cm x 78.7 cm)

### Accessories

- › 251 – Bierhoff Knee Crutches (pair)
- › 582 – Hand Control
- › 585 – Armrests
- › 590 – Welch Allyn Light Bracket  
capabilities for #48955 (patient left)
- › 591 – Welch Allyn Light Bracket  
capabilities for #48955 (patient right)



Model 3503-500-300  
includes Classic upholstery

### Upholstery and Colors

Soft Linen	Smoky Cashmere	Warm Sand	Creamy Latte
True Graphite	Classic Black	Fine Wine	Blue Skies
Chocolate Truffle	Lemon Meringue	Cherry Blossom	Morning Fog
Lakeside Blue	Twilight Blue	Mint Leaf	Deep Forest



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

©2023 United Metal Fabricators Inc. dba UMF Medical 07.23 All Rights Reserved. PreFixx® is a registered trademark of Omnova Solutions. FusionONE is a trademark of United Metal Fabricators. 5-1372LT Rev A.



Ph: +1 814.266.8726 | Fax: +1 814.266.1870  
1316 Eisenhower Blvd., Johnstown, PA 15904 USA

©2023 United Metal Fabricators Inc. dba UMF Medical

**UMFmedical.com**