

Four Hook, Stainless Steel I.V. Pole Mounted on Steel Base

Model SS8344

The UMF Medical SS844 I.V. Pole has 4 hooks and a 22.5" wheel base diameter. The unit is mounted on 4 legged 3/8" steel base equipped with 3" swivel casters. The telescopic stainless steel rod is secured in place by a positive-locking hand wheel. The I.V. pole comes with a 3-year limited warranty.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Specifications

- › 4 Hooks
- › Wheel base diameter 22.5" (57.2cm)
- › Caster Size 3" (7.6cm)
- › 4 leg construction using 3/8" steel
- › Friction knob to secure height
- › Compressed length 68" (172.7cm)
- › Extended length 96" (243.8cm)
- › Outer base tube is 63.25" (160.7cm) long
- › Outer base tube is 1" diameter stainless steel
- › Maximum weight capacity is 40lb (18kg) when extended to 6ft height
- › 3 year limited warranty