

power200+

power exam table

Power Hi-Lo, Power Back Exam Table with OneTouch Patient Positioning® System, Hand and Foot Control

Model 4040-650-200

For examinations

The power200+ is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers.

Designed for exams, the power200+ meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.



PROTECTS > patients from infection > facilities from damage > staff from injury

Key features at a glance

- > Power hi-lo and power backrest
- > Two function foot control with hi-lo and power back
- > Hand Control with OneTouch Patient Positioning® 2 programmable positions & return-to-home
- > Dual electro-mechanical actuator system for efficient speed
- > 18"-21" paper roll holder (53.3cm)
- > 650lb (294.84kg) patient weight capacity
- > Low-access ADA height 18.5" (±0.5") - 46.9cm (±1.27cm)
- > Pelvic tilt (7°)
- > Drain pan: stainless steel 1.5" deep (3.81cm)
- > Front storage drawer
- > Side storage drawer – reversible
- > Removable, seamless upholstery with Prefix® Protective Finish
- > Stirrups with four lateral positions, adjustable length
- > Front leg cushion for patient comfort
- > Heavy gauge steel extension leg pull-out to accommodate all patient heights
- > Hospital grade receptacle 110V with child safety cover (patient left)
- > 3 year limited warranty
- > Durable galvanized steel construction
- > Easy-clean powder coat enamel surfaces

Wireless Hand & Foot Controls Available



Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm) at seat max. 37" at seat (93.9cm)
- › Width: 28" (71.12cm)
- › Length: 57" – 76" with foot section extended (144.7cm - 193cm)
- › Leg rest: 19" W x 13" D (48.26cm W x 33cm D)
- › Front storage drawer: 19" W x 12" D (48.26cm W x 30.48cm D)
- › Stirrups: four lateral positions, adjustable length
- › Dual electro-mechanical actuator system for efficient speed
- › Electrical: 120 VAC, 60 HZ
- › International model: 4040-650-201 220/240V, 50/60Hz
- › Net weight: 445lb (202kg)
- › Shipping weight: 485lb (220kg)
- › Shipping dimensions: L64 x W30 x H31 (163cm x 76cm x 79cm)

Accessories

- › 042 – Patient assist arms complete (pair) (field installable)
- › 045 – Universal Accessory Rail-Seat Section
- › 234 – Two Hook IV pole kit (includes IV pole and clamp)
- › 235 – Four Hook IV pole kit (includes IV pole and clamp)
- › 236 – Fixed armboard package (EA)
- › 237 – Articulating armboard package (EA)
- › 251 – Bierhoff knee crutches (pair)
- › 442 – Two-function wireless foot control
- › 443 – Two-function wireless hand control
- › 467 – Side Mounting plate for Welch Allyn Light Bracket #48955 (Patient Right), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 468 – Side Mounting plate for Welch Allyn Light Bracket #48955 (PatientLeft), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 471 – Front Mounting Plate for Welch Allyn Light Bracket #48955 (Patient Left), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 472 – Front Mounting Plate for Welch Allyn Light Bracket #48955 (Patient Right), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 822 – Removable stainless steel drain pan w/spout
- › 870 – 110v drawer warmer (factory installed)
- › 872 – 220v drawer warmer (factory installed) (international)
- › 873 – 220v outlet package (factory installed) (international)
- › 910 – Easy power table installation kit

Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA