

power4010^P

head-centric procedure chair

Procedure Chair with OneTouch Patient Positioning[®] System

For head-centric practices

The power4010p is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers, such as ENT, Dermatology, and Plastic Surgery professionals.

Designed for head-centric procedures, the power4010p meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.



Model 4010-650-200

PROTECTS > patients from infection > facilities from damage > staff from injury

Key features at a glance

- > Chair designed for head-centric practices
- > 650 lb. (294.84 kg) patient weight capacity
- > Low Access ADA Height 18.5" (±0.5") 46.9cm (±1.27cm)
- > 4 Power Functions: hi-lo, back, tilt, leg lift & extension
- > Hand control with OneTouch Patient Positioning[®] System, 2 programmable positions & return-to-home
- > Foot control with OneTouch Patient Positioning[®] System, 2 programmable positions and return-to-home
- > Double-articulating headrest with adjustable height and interchangeable headrest options
- > Slim back cushion for accessibility
- > Removable, seamless upholstery with Prefix[®] Protective Finish
- > Trendelenburg, Lithotomy, Semi-Fowler positions
- > Hospital-Grade electrical receptacle with child safety cover (patient left)
- > Concealed 18" paper roll holder
- > 3 year limited warranty

Hand Control



Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm) at seat max. 37" at seat 93.9cm
- › Seat width: 28" (71.12cm)
- › Back width: 24.5" (62.23cm)
- › Length: 70" – 80" with foot section extended (177.8cm – 203.2cm)
- › Leg rest: 28" W x 13" D (71.12cm W x 33.02cm D)
- › Dual electro-mechanical actuator system for efficient speed
- › Electrical: 120 VAC, 60HZ
- › International model: 4010-650-101: 220/240V, 50/60Hz
- › Net weight: 410lb (185kg)
- › Shipping weight: 450lb (204kg)
- › Shipping dimensions: L64 x W30 x H41 (163cm x 76cm x 104cm)

Accessories

- › 047 – Adjustable / self-leveling arms (pair)
- › 050 – Removable U-shaped headrest
- › 051 – Removable rectangular headrest
- › 234 – 2 Hook IV Pole Kit (includes clamp)
- › 235 – 4 Hook IV Pole Kit (includes clamp)
- › 238 – Fixed armboard package (EA)
- › 239 – Articulating armboard package (EA)
- › 473 – Welch Allyn light bkt capabilities for #48955 (patient left)
- › 474 – Welch Allyn light bkt capabilities for #48955 (patient right)
- › 830 – Patient safety strap (factory installed only)
- › 910 – Easy power table installation kit



Shown with options (047 Self-Leveling Arms)

Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2020 United Metal Fabricators Inc. dba UMF Medical 05.23. All Rights Reserved. OneTouch Patient Positioning® is a registered trademark of United Metal Fabricators Inc. 5-1326LT Rev A