

# proglide4010


head centric procedure chair

## Procedure Chair with OneTouch WheelBase® System, OneTouch Patient Positioning® System, Hand and Foot Control

### For head-centric practices

The proglide4010 is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers, such as ENT, Dermatology, and Plastic Surgery professionals.

Designed for head-centric procedures, the proglide4010 meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.



Model 4010-650-300

**PROTECTS** > patients from infection > facilities from damage > staff from injury

### Key features at a glance

- > Chair designed for head-centric practices
- > OneTouch Wheelbase®
- > 650 lb. (294.84 kg) patient weight capacity
- > Low access ADA height 18.5" (±0.5") 46.9cm (±1.27cm)
- > 4 Power functions: hi-lo, back, tilt, leg lift & extension
- > Hand control with OneTouch
- > Patient Positioning® System, 2 programmable positions & return-to-home
- > Foot control with OneTouch Patient Positioning® System, 2 programmable positions and return-to-home
- > Double-articulating headrest with height adjustability and interchangeable headrest options
- > Slim back cushion for accessibility
- > Removable seamless upholstered cushion
- > Trendelenberg, Lithotomy & Semi-Fowler positions
- > Hospital-Grade electrical receptacle with child safety cover (patient left)
- > Concealed 18" paper roll holder
- > 3 year limited warranty

Hand Control



## Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm)  
at seat max. 37" at seat (93.9cm)
- › Seat width: 28" (71.12cm)
- › Back width: 24.5" (62.23cm)
- › Length: 70" – 80" with foot section  
extended (177.8cm – 203.2cm)
- › Leg rest: 28" W x 13" D (71.12cm  
W x 33.02cm D)
- › Dual electro-mechanical actuator  
system for efficient speed
- › Electrical: 120 VAC, 60HZ , 5 AMP
- › International model: 4010-650-301:  
220/240V, 50/60Hz
- › Net weight: 450lb (204kg)
- › Shipping weight: 490lb (222kg)
- › Shipping dimensions: L64 x W30 x H41  
(163cm x 76cm x 104cm)

## Accessories

- › 047 – Adjustable / self-leveling arms (pair)
- › 050 – Removable U-shaped headrest
- › 051 – Removable rectangular headrest
- › 238 – Fixed armboard package (EA)
- › 239 – Articulating armboard package (EA)
- › 473 – Welch Allyn light bkt capabilities  
for #48955 (patient left)
- › 474 – Welch Allyn light bkt capabilities  
for #48955 (patient right)
- › 830 – Patient safety strap  
(factory installed only)
- › 910 – Easy power table installation kit
- › Fire retardant foam (cal133) available



Shown with  
options (047  
Self-Leveling Arms)

## Colors



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

©2020 United Metal Fabricators Inc. dba UMF Medical 06.19 All Rights Reserved. OneTouch Patient Positioning® is a registered trademark of United Metal Fabricators Inc. 5-1327LT Rev A