

# proglide300<sup>+</sup>

power exam table

## Power Hi-Lo, Power Back Exam Chair with OneTouch WheelBase<sup>®</sup> System, OneTouch Patient Positioning System, Hand and Foot Control

Model 4040-650-200



### For examinations

The proglide300+ is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers. Designed for exams, the proglide300+ meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.

**PROTECTS** > patients from infection > facilities from damage > staff from injury

### Key features at a glance

- > OneTouch WheelBase<sup>®</sup> integrated central-locking caster system
- > Hand control with OneTouch Patient Positioning<sup>®</sup>, 2 programmable positions & return-to-home
- > Two function foot control with hi-lo and power back
- > Dual electro-mechanical actuator system for efficient speed
- > 18"-21" paper roll holder (53.3cm)
- > 650lb (294.84kg) patient weight capacity
- > Low-access ADA height 18.5" (±0.5") - 46.9cm (±1.27cm)
- > Pelvic tilt (7°)
- > Drain pan: stainless steel 1.5" deep (3.81cm)
- > Front storage drawer
- > Side storage drawer – reversible
- > Removable, seamless upholstered top
- > Stirrups with four lateral positions, adjustable length
- > Front leg cushion for patient comfort
- > Heavy gauge steel extension leg pull-out to accommodate all patient heights
- > Drawer warmer
- > Hospital grade receptacle 110V with child safety cover (patient left)
- > 3 year limited warranty
- > Durable galvanized steel construction
- > Easy-clean powder coat enamel surfaces

## Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm) at seat max. 37" at seat (93.9cm)
- › Width: 28" (71.12cm)
- › Length: 57" – 76" with foot section extended (144.7cm - 193cm)
- › Leg rest: 19" W x 13" D (48.26cm W x 33cm D)
- › Front storage drawer: 19" W x 12" D (48.26cm W x 30.48cm D)
- › Stirrups: four lateral positions, adjustable length
- › Dual electro-mechanical actuator system for efficient speed
- › Electrical: 120 VAC, 60 HZ
- › International model 4040-650-301: 220/240V, 50/60Hz
- › Net weight: 445lb (202kg)
- › Shipping weight: 485lb (220kg)
- › Shipping dimensions: L64 x W30 x H31 (163cm x 76cm x 79cm)

## Accessories

- › 022 – Chair-to-bed elevated leg pad creates flat transfer surface to meet US access board standards
- › 042 – Patient assist arms complete (pair) (field installable)
- › 236 – Fixed armboard package (EA)
- › 237 – Articulating armboard package (EA)
- › 251 – Bierhoff knee crutches (pair)
- › 471 – Welch Allyn light bkt capabilities for #48955 (patient left)
- › 472 – Welch Allyn light bkt capabilities for #48955 (patient right)
- › 822 – Removable stainless steel drain pan w/spout
- › 872 – 220v drawer warmer (factory installed) (international)
- › 873 – 220v outlet package (factory installed) (international)
- › 910 – Easy power table installation kit
- › Fire retardant foam (cal133) available

## Colors



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

©2019 United Metal Fabricators Inc. dba UMF Medical 06.19 All Rights Reserved. OneTouch Patient Positioning® is a registered trademark of United Metal Fabricators Inc. 5-1333LT Rev A