

proglide300

power exam table

Power Hi-Lo, Manual Back Exam Table with OneTouch WheelBase[®], OneTouch Patient Positioning[®] System, Hand and Foot Control

Model 4070-650-300



For examinations

The proglide300 is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers.

PROTECTS > patients from infection > facilities from damage > staff from injury

Key features at a glance

- > OneTouch WheelBase[®] integrated central-locking caster system
- > Hand control with OneTouch Patient Positioning[®], 2 programmable positions & return-to-home
- > One function foot control
- > Power hi-lo with smooth moving articulating back adjustment
- > Dual electro-mechanical actuator system for efficient speed
- > 18"-21" paper roll holder (53.3cm)
- > 650lb (294.84kg) patient weight capacity
- > Seat height: 18.5" (±0.5")-46.9cm (±1.27cm)
- > Pelvic tilt (7°)
- > Drain pan: stainless steel 1.5" deep (3.81cm)
- > Front storage drawer
- > Side storage drawer – reversible
- > Removable, seamless upholstery with Prefixx[®] Protective Finish
- > Stirrups with four lateral positions, adjustable length
- > Front leg cushion for patient comfort
- > Heavy gauge steel extension leg pull-out to accommodate all patient heights
- > Drawer warmer
- > Hospital grade receptacle 110V with child safety cover (patient left)
- > 3 year limited warranty
- > Durable galvanized steel construction
- > Easy-clean powder coat enamel surfaces

Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm) at seat max. 37" at seat (93.9cm)
- › Width: 28" (71.12cm)
- › Length: 57" - 76" with foot section extended (144.7cm - 193cm)
- › Leg rest: 19" W x 13" D (48.26cm W x 33cm D)
- › Front storage drawer: 19" W x 12" D (48.26cm W x 30.48cm D)
- › Stirrups: four lateral positions, adjustable length
- › Dual electro-mechanical actuator system for efficient speed
- › Electrical: 120 VAC, 60 HZ
- › International model 4070-650-301: 220/240V, 50/60Hz
- › Net weight: 470lb (213kg)
- › Shipping weight: 510lb (231kg)
- › Shipping dimensions: L64xW30xH31 (163cm x 76cm x 79cm)

Accessories

- › 042 - Patient assist arms complete (pair) field installable
- › 045 - Universal Accessory Rail-Seat Section
- › 234 - Two Hook IV pole kit (includes IV pole and clamp)
- › 235 - Four Hook IV pole kit (includes IV pole and clamp)
- › 236 - Fixed armboard package (EA)
- › 237 - Articulating armboard package (EA)
- › 251 - Bierhoff knee crutches (pair)
- › 467 - Side Mounting plate for Welch Allyn Light Bracket #48955 (Patient Right), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 468 - Side Mounting plate for Welch Allyn Light Bracket #48955 (PatientLeft), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 471 - Front Mounting Plate for Welch Allyn Light Bracket #48955 (Patient Left), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 472 - Front Mounting Plate for Welch Allyn Light Bracket #48955 (Patient Right), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 822 - Removable stainless steel drain pan W/spout
- › 872 - 220v drawer warmer (factory installed) (international)
- › 873 - 220v Outlet Package (factory installed) (international)
- › 910 - Easy power table installation kit

Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2026 United Metal Fabricators Inc. dba UMF Medical 06.26 All Rights Reserved. OneTouch Patient Positioning® is a registered trademark of United Metal Fabricators Inc. 5-1336LT Rev B