

proglide300

power exam table

Power Hi-Lo, Manual Back Exam Table with OneTouch WheelBase[®], OneTouch Patient Positioning[®] System, Hand and Foot Control

Model 4070-650-300



For examinations

The proglide300 is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers. Designed for exams, the proglide300 meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.

PROTECTS > patients from infection > facilities from damage > staff from injury

Key features at a glance

- > OneTouch WheelBase[®] integrated central-locking caster system
- > Hand control with OneTouch Patient Positioning[®], 2 programmable positions & return-to-home
- > One function foot control
- > Power hi-lo with smooth moving articulating back adjustment
- > Dual electro-mechanical actuator system for efficient speed
- > 18"-21" paper roll holder (53.3cm)
- > 650lb (294.84kg) patient weight capacity
- > Low-access ADA height 18.5" (±0.5") - 46.9cm (±1.27cm)
- > Pelvic tilt (7°)
- > Drain pan: stainless steel 1.5" deep (3.81cm)
- > Front storage drawer
- > Side storage drawer – reversible
- > Removable, seamless upholstery with Prefix[®] Protective Finish
- > Stirrups with four lateral positions, adjustable length
- > Front leg cushion for patient comfort
- > Heavy gauge steel extension leg pull-out to accommodate all patient heights
- > Drawer warmer
- > Hospital grade receptacle 110V with child safety cover (patient left)
- > 3 year limited warranty
- > Durable galvanized steel construction
- > Easy-clean powder coat enamel surfaces

Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm) at seat max. 37" at seat (93.9cm)
- › Width: 28" (71.12cm)
- › Length: 57" - 76" with foot section extended (144.7cm - 193cm)
- › Leg rest: 19" W x 13" D (48.26cm W x 33cm D)
- › Front storage drawer: 19" W x 12" D (48.26cm W x 30.48cm D)
- › Stirrups: four lateral positions, adjustable length
- › Dual electro-mechanical actuator system for efficient speed
- › Electrical: 120 VAC, 60 HZ
- › International model 4070-650-301: 220/240V, 50/60Hz
- › Net weight: 470lb (213kg)
- › Shipping weight: 510lb (231kg)
- › Shipping dimensions: L64xW30xH31 (163cm x 76cm x 79cm)

Accessories

- › 042 - Patient assist arms complete (pair) field installable
- › 045 - Universal Accessory Rail-Seat Section
- › 234 - Two Hook IV pole kit (includes IV pole and clamp)
- › 235 - Four Hook IV pole kit (includes IV pole and clamp)
- › 236 - Fixed armboard package (EA)
- › 237 - Articulating armboard package (EA)
- › 251 - Bierhoff knee crutches (pair)
- › 467 - Side Mounting plate for Welch Allyn Light Bracket #48955 (Patient Right), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 468 - Side Mounting plate for Welch Allyn Light Bracket #48955 (PatientLeft), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 471 - Front Mounting Plate for Welch Allyn Light Bracket #48955 (Patient Left), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 472 - Front Mounting Plate for Welch Allyn Light Bracket #48955 (Patient Right), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 822 - Removable stainless steel drain pan W/spout
- › 872 - 220v drawer warmer (factory installed) (international)
- › 873 - 220v Outlet Package (factory installed) (international)
- › 910 - Easy power table installation kit

Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2020 United Metal Fabricators Inc. dba UMF Medical 05.23 All Rights Reserved. OneTouch Patient Positioning® is a registered trademark of United Metal Fabricators Inc. 5-1336LT Rev A