

# power100

power exam table

## Power Hi-Low, Manual Back Exam Table with Foot Control

### For examinations

The power100 is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers.

Designed for exams, the power100 meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.



Model 4070-650-100

**PROTECTS** > patients from infection > facilities from damage > staff from injury

### Key features at a glance

- > Power Hi-Lo with smooth moving articulating back adjustment
- > Dual electro-mechanical actuator system for efficient speed
- > One function foot control
- > 18"-21" Paper roll holder (53.3cm)
- > 650lb (294.84kg) patient weight capacity
- > Low-access ADA height 18.5" ( $\pm 0.5$ ") - 46.9cm ( $\pm 1.27$ cm)
- > Pelvic tilt (7°)
- > Drain pan: stainless steel 1.5" deep (3.81cm)
- > Front storage drawer
- > Removable, seamless upholstered cushion
- > Stirrups with four lateral positions, adjustable length
- > Front leg cushion for patient comfort
- > Heavy gauge steel extension leg pull-out to accommodate all patient heights
- > 3 year limited warranty
- > Durable galvanized steel construction
- > Easy-clean powder coat enamel surfaces

## Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm) at seat max. 37" at seat (93.9cm)
- › Width: 28" (71.12cm)
- › Length: 57" – 76" with foot section extended (144.7cm - 193cm)
- › Leg Rest: 19" W x 13" D (48.26cm W x 33cm D)
- › Front Storage Drawer: 19" W x 12" D (48.26cm W x 30.48cm D)
- › Stirrups: four lateral positions, adjustable length
- › Dual electro-mechanical actuator system for efficient speed
- › Electrical: 120 VAC, 60 HZ
- › International model: 4070-650-101 220/240V, 50/60Hz
- › Net weight: 400lb (181kg)
- › Shipping weight: 445lb (202kg)
- › Shipping dimensions: L64 x W30 x H31 (163cm x 76cm x 79cm)

## Accessories

- › 022 – Chair-to-bed elevated leg pad creates flat transfer surface to meet US access board standards
- › 042 – Patient assist arms complete (pair) field installable
- › 236 – Fixed armboard package (EA)
- › 237 – Articulating armboard package (EA)
- › 251 – Bierhoff knee crutches (pair)
- › 471 – Welch Allyn light bkt capabilities for #48955 (patient left)
- › 472 – Welch Allyn light bkt capabilities for #48955 (patient right)
- › 822 – Removable stainless steel drain pan W/spout
- › 870 – 110v drawer warmer (factory installed)
- › 871 – 110v outlet package (factory installed)
- › 872 – 220v drawer warmer (factory installed) (international)
- › 873 – 220v outlet package (factory installed) (international)
- › 910 – Easy power table installation kit
- › Fire retardant foam (cal133) available

## Colors



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

©2020 United Metal Fabricators Inc. dba UMF Medical 06.19 All Rights Reserved. OneTouch Patient Positioning® is a registered trademark of United Metal Fabricators Inc. 5-1334LT Rev A