

# power100

power exam table

## Power Hi-Low, Manual Back Exam Table with Foot Control

### For examinations

The power100 is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers.



Model 4070-650-100

**PROTECTS** > patients from infection > facilities from damage > staff from injury

### Key features at a glance

- > Power Hi-Lo with smooth moving articulating back adjustment
- > Dual electro-mechanical actuator system for efficient speed
- > One function foot control
- > 18"-21" Paper roll holder (53.3cm)
- > 650lb (294.84kg) patient weight capacity
- > Seat height: 18.5" (±0.5") - 46.9cm (±1.27cm)
- > Pelvic tilt (7°)
- > Drain pan: stainless steel 1.5" deep (3.81cm)
- > Front storage drawer
- > Removable, seamless upholstery with Prefixx<sup>®</sup> Protective Finish
- > Stirrups with four lateral positions, adjustable length
- > Front leg cushion for patient comfort
- > Heavy gauge steel extension leg pull-out to accommodate all patient heights
- > 3 year limited warranty
- > Durable galvanized steel construction
- > Easy-clean powder coat enamel surfaces

## Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm) at seat max. 37" at seat (93.9cm)
- › Width: 28" (71.12cm)
- › Length: 57" – 76" with foot section extended (144.7cm - 193cm)
- › Leg Rest: 19" W x 13" D (48.26cm W x 33cm D)
- › Front Storage Drawer: 19" W x 12" D (48.26cm W x 30.48cm D)
- › Stirrups: four lateral positions, adjustable length
- › Dual electro-mechanical actuator system for efficient speed
- › Electrical: 120 VAC, 60 HZ
- › International model: 4070-650-101 220/240V, 50/60Hz
- › Net weight: 400lb (181kg)
- › Shipping weight: 445lb (202kg)
- › Shipping dimensions: L64 x W30 x H31 (163cm x 76cm x 79cm)

## Accessories

- › 042 – Patient assist arms complete (pair) field installable
- › 045 – Universal Accessory Rail-Seat Section
- › 236 – Fixed armboard package (EA)
- › 237 – Articulating armboard package (EA)
- › 251 – Bierhoff knee crutches (pair)
- › 467 – Side Mounting plate for Welch Allyn Light Bracket #48955 (Patient Right), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 468 – Side Mounting plate for Welch Allyn Light Bracket #48955 (PatientLeft), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 471 – Front Mounting Plate for Welch Allyn Light Bracket #48955 (Patient Left), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 472 – Front Mounting Plate for Welch Allyn Light Bracket #48955 (Patient Right), Compatible with WA Green Series Light 48810, Field Installable, for Series 4040-650, 4070-650
- › 822 – Removable stainless steel drain pan W/spout
- › 870 – 110v drawer warmer (factory installed)
- › 871 – 110v outlet package (factory installed)
- › 872 – 220v drawer warmer (factory installed) (international)
- › 873 – 220v outlet package (factory installed) (international)
- › 910 – Easy power table installation kit

## Colors



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

©2026 United Metal Fabricators Inc. dba UMF Medical 06.26 All Rights Reserved. OneTouch Patient Positioning® is a registered trademark of United Metal Fabricators Inc. 5-1334LT Rev B