

power4011^P

ultra procedure chair

Ultra Procedure Chair with OneTouch Patient Positioning[®] System Hand and Foot Control

For procedures

The power4011p is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers, such as OB/GYN and Oncology professionals.

Designed for procedures, the power4011p meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.



Model 4011-650-200

PROTECTS > patients from infection > facilities from damage > staff from injury

Key features at a glance

- > 650 lb. (294.84kg) patient weight capacity
- > Low access ADA height 18.5" (±0.5") 46.9cm (±1.27cm)
- > 4 power functions: hi-lo, back, tilt, leg lift & extension
- > Hand control with OneTouch Patient Positioning[®] System, 2 programmable positions & return-to-home
- > Foot control with OneTouch Patient Positioning[®] System, 2 programmable positions & return-to-home
- > Removable seamless upholstered top
- > Double-articulating headrest with adjustable height
- > Interchangeable headrest options
- > Stirrups
- > Stainless steel drain pan
- > Trendelenburg, Lithotomy, Semi-Fowler positions
- > Hospital-grade electrical receptacle with child safety cover (patient left)
- > Concealed paper roll holder
- > 3 year limited warranty



Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm)
at seat max. 37" at seat (93.9cm)
- › Seat width: 28" (71.12cm)
- › Back width: 26.5" (67.31cm)
- › Length: 73"-80" with foot section
extended (185.42cm – 203.2cm)
- › Leg rest: 28" W x 11.5" D
(71.12cm W x 29.21cm D)
- › Stirrups: four lateral positions,
adjustable length
- › Dual electro-mechanical actuator
system for efficient speed
- › Electrical: 120 VAC, 60HZ
- › International model: 4011-650-201:
220/240V, 50/60Hz
- › Net weight: 465lb (211kg)
- › Shipping weight: 505lb (229kg)
- › Shipping dimensions: L64 x W30 x H41
(163cm x 76cm x 104cm)

Accessories

- › 044 – Adjustable / self-leveling arms (pair)
- › 050 – Removable U-shaped headrest
- › 051 – Removable rectangular headrest
- › 238 – Fixed armboard package (EA)
- › 239 – Articulating armboard package (EA)
- › 473 – Welch Allyn light bkt capabilities
for #48955 (patient left)
- › 474 – Welch Allyn light bkt capabilities
for #48955 (patient right)
- › 830 – Patient safety strap
(factory installed only)
- › 873 – 220v outlet package
(factory installed) (international)
- › 874 – 220v power cord (international)
- › 910 – Easy power table installation kit
- › Fire retardant foam (cal133) available



Shown with
options (044
Self-Leveling Arms)

Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2020 United Metal Fabricators Inc. dba UMF Medical 06.19 All Rights Reserved. OneTouch Patient Positioning® is a registered trademark of United Metal Fabricators Inc. 5-1329LT Rev A