

# proglide4011

ultra procedure chair

Model 4011-650-300  
Shown with optional  
Self-Leveling Arms

## Ultra Procedure Chair with OneTouch WheelBase<sup>®</sup>, OneTouch Patient Positioning<sup>®</sup> System Hand and Foot Control

### For procedures

The proglide4011 is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers, such as OB/GYN and Oncology professionals. Designed for, the Proglide4011 accommodates the needs of most patients.



**PROTECTS** > patients from infection > facilities from damage > staff from injury

### Key features at a glance

- > OneTouch WheelBase<sup>®</sup> integrated central-locking caster system
- > 650 lb. (294.84kg) patient weight capacity
- > Low height 18.5" (±0.5") 46.9cm (±1.27cm)
- > 4 power functions: hi-lo, back, tilt, leg lift & extension
- > Hand control with OneTouch Patient Positioning<sup>®</sup> System, 2 programmable positions & return-to-home
- > Foot control with OneTouch Patient Positioning<sup>®</sup> System, 2 programmable positions & return-to-home
- > Removable, seamless upholstered top with Prefix<sup>®</sup> Protective Finish
- > Double-articulating headrest with adjustable height
- > Stirrups
- > Interchangeable headrest options
- > Stainless steel drain pan
- > Trendelenburg, Lithotomy, Semi-Fowler positions
- > Hospital-grade electrical receptacle with child safety cover (patient left)
- > Concealed paper roll holder
- > 3 year limited warranty

Hand Control



## Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm)  
at seat max. 37" at seat (93.9cm)
- › Seat width: 28" (71.12cm)
- › Back width: 26.5" (67.31cm)
- › Length: 73"-80" with foot section  
extended (185.42cm – 203.2cm)
- › Leg rest: 28" W x 11.5" D  
(71.12cm W x 29.21cm D)
- › Stirrups: four lateral  
positions, adjustable length
- › Dual electro-mechanical actuator  
system for efficient speed
- › Electrical: 120 VAC, 60HZ
- › International model: 4011-650-301:  
220/240V, 50/60Hz
- › Net weight: 495lb (224kg)
- › Shipping weight: 535lb (242kg)
- › Shipping dimensions: L64 x W30 x H41  
(163cm x 76cm x 104cm)

## Accessories

- › 044 – Adjustable / self-leveling arms (pair)
- › 050 – Removable U-shaped headrest
- › 051 – Removable rectangular headrest
- › 234 – 2 Hook IV Pole Kit (includes clamp)
- › 235 – 4 Hook IV Pole Kit (includes clamp)
- › 238 – Fixed armboard package (EA)
- › 239 – Articulating armboard package (EA)
- › 251 – Bierhoff knee crutches (pair)
- › 473 – Welch Allyn light bkt capabilities  
for #48955 (patient left)
- › 474 – Welch Allyn light bkt capabilities  
for #48955 (patient right)
- › 830 – Patient safety strap  
(factory installed only)
- › 910 – Easy power table installation kit



Shown with  
options (044  
Self-Leveling Arms)

## Colors



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

©2020 United Metal Fabricators Inc. dba UMF Medical 05.23. All Rights Reserved. OneTouch Patient Positioning® is a registered trademark of United Metal Fabricators Inc. 5-1330LT Rev A