

power5016^P

podiatry/wound care procedure chair

Model 5016-650-200
Shown with optional
Self-Leveling Arms

Procedure Chair, Power Leg Lift, Tilt, Hi-Lo & Back, OneTouch Patient Positioning[®] System Hand & Foot Controls

For procedures

The power5016p is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers. Designed for podiatry and wound care procedures, the Power 5016p accommodates the needs of most patients.



PROTECTS

- ▶ patients from infection
- ▶ facilities from damage
- ▶ staff from injury

Key features at a glance

- › OneTouch Patient Positioning[®] programmable foot control for efficiencies
- › OneTouch Patient Positioning[®] programmable hand control
- › 4 Power Functions: hi-lo, back, tilt, leg lift & extension
- › Power leg lift makes it easier for patients, easier for staff
- › Front access to chair
- › Removable, Seamless Upholstered Top with Prefixx[®] Protective Finish
- › Adjustable articulating head rest provides ultimate comfort for patients
- › 650 lb. patient weight capacity
- › Low height 18.5" (±0.5") 46.9cm (±1.27cm)
- › Hospital-grade electrical receptacle 110v with child safety cover (patient left)
- › Concealed 18" paper roll holder
- › 3 year limited warranty
- › Interchangeable headrest options

Hand
Control



Specifications

- › Height: min. 18.5" (± .5") 46.9cm (±1.27cm) at seat max. 37" at seat (93.9cm)
- › Seat width: 28" (71.12cm)
- › Back width: 26.5" (67.31cm)
- › Length: 73-80" with foot section extended (185.42cm – 203.2cm)
- › Leg rest: 28" W x 13" D (71.12cm W x 33.02cm D)
- › Dual electro-mechanical actuator system for efficient speed
- › Electrical: 120 VAC, 60HZ - 220/240V, 50/60Hz
- › International model (5016-650-201)
- › Net weight: 421lb (191kg)
- › Shipping weight: 461lb (209kg)
- › Shipping dimensions: L64xW30xH41 (163cm x 76cm x 104cm)

Accessories

- › 044 – Adjustable / self-leveling arms (pair)
- › 045 – Universal Accessory Rail Installed on Seat Section
- › 050 – Removable U-shaped headrest
- › 051 – Removable rectangular headrest
- › 234 – 2 Hook IV Pole Kit (includes clamp)
- › 235 – 4 Hook IV Pole Kit (includes clamp)
- › 238 – Fixed armboard package (EA)
- › 239 – Articulating armboard package (EA)
- › 473 – Welch Allyn light bkt capabilities for #48955 (patient left)
- › 474 – Welch Allyn light bkt capabilities for #48955 (patient right)
- › 830 – Patient safety strap (factory installed only)
- › 873 – 220v outlet package (factory installed) (international)
- › 910 – Easy power table installation kit

Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

©2024 United Metal Fabricators Inc. dba UMF Medical 12.24. All Rights Reserved. OneTouch Patient Positioning® is a registered trademark of United Metal Fabricators Inc. 5-1338LT Rev A