

# proglide5016

podiatry/wound care procedure chair

## Procedure Chair, Power Leg Lift, Tilt, Hi-Lo & Back, OneTouch Patient Positioning® System Hand & Foot Controls

### For procedures

The proglide5016 is a product of extensive research and clinical feedback, developed to improve patient care and streamline the work of healthcare providers.

Designed for podiatry and wound care procedures, the proglide5016 meets ADA and U.S. Access Board regulations and accommodates the needs of most patients.



Model 5016-650-300

**PROTECTS** > patients from infection > facilities from damage > staff from injury

### Key features at a glance

- > OneTouch WheelBase® integrated central-locking caster system
- > One Touch Patient Positioning® System foot control for efficiencies
- > OneTouch Patient Positioning® System hand control
- > 4 Power functions: hi-Lo, back, tilt, leg lift & extension
- > Power leg lift makes it easier for patients, easier for staff
- > Front access to chair
- > Ultra-comfort seat and back cushion
- > Adjustable articulating head rest provides ultimate comfort for patients
- > 650 lb. patient weight capacity
- > Low access ADA height 18.5" ( ±0.5" ) 46.9cm ( ±1.27cm)
- > Hospital-grade electrical receptacle 110v with child safety cover (patient left)
- > Concealed 18" paper roll holder
- > 3 year limited warranty
- > Interchangeable headrest options



## Specifications

- › Height: Min. 18.5" (± .5") 46.9cm (±1.27cm) at seat. Max. 37" at seat (93.9cm)
- › Seat width: 28" (71.12cm)
- › Back width: 26.5" (67.31cm)
- › Length: 73-80" with foot section extended (185.42cm - 203.2cm)
- › Leg rest: 28" W x 13" D (71.12cm W x 33.02cm D)
- › Dual electro-mechanical actuator system for efficient speed
- › Electrical: 120 VAC, 60HZ - 220/240V, 50/60Hz
- › International model (5016-650-301)
- › Net weight: 421lb (191kg)
- › Shipping weight: 461lb (209kg)
- › Shipping dimensions: L64xW30xH41 (163cm x 76cm x 104cm)

## Accessories

- › 044 - Adjustable / self-leveling arms (pair)
- › 050 - Removable U-shaped headrest
- › 051 - Removable rectangular headrest
- › 238 - Fixed armboard package (EA)
- › 239 - Articulating armboard package (EA)
- › 235 - Hand surgery armboard (EA)
- › 473 - Welch Allyn light bkt capabilities for #48955 (patient left)
- › 474 - Welch Allyn light bkt capabilities for #48955 (patient right)
- › 830 - Patient safety strap (factory installed only)
- › 873 - 220v outlet package (factory installed) (international)
- › 910 - Easy power table installation kit
- › Fire retardant foam (cal133) available

## Colors



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

©2020 United Metal Fabricators Inc. dba UMF Medical 06.19 All Rights Reserved. OneTouch Patient Positioning® is a registered trademark of United Metal Fabricators Inc. 5-1339LT Rev A