



## Portable Emergency Spring Bed Frame

### Model 8560

UMF Medical Portable Emergency Spring Bed Frames have a comfort-spring frame and require minimal, tool-free assembly and disassembly for quick setup and teardown. Available with or without mattress, UMF Medical Emergency Beds are ideal for fast response, temporary stays, and can be shipped in large quantities. Spring bed frame ships assembled, two leg sections ship separately. Mattress sold separately. Easy to transport disassembled.



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

### Key features at a glance

- › Comfort-spring, platform spring bed frame
- › Minimal, tool-free assembly and disassembly for quick, efficient setup and teardown in emergent situations
- › Lightweight and easy to carry disassembled.
- › The durable, powder-coated steel frame provides maximum support
- › 17-inch frame height allows easy access for a varied patient population
- › 15 inches of under bed storage for medical equipment and supplies.
- › Spring frame ships assembled, two leg sections ship separately.
- › Mattress sold separately.

### Specifications

Weight capacity: 350 lbs.  
 Dimensions: 30" w x 74" L x 17" H  
 Weight: 35lbs (flat surface) 10lbs (legs)  
 Packaging Quantities: 20 beds/pallet

### Accessories

- › 8570 – 4in. mattress