



## Stainless Steel Instrument Table with Shelf

### Model SS8000

UMF Medical stainless steel instrument tables are made with 20 gauge stainless steel shelves and 16 gauge stainless steel legs. All legs come with swivel casters and a 3-year limited warranty.



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

#### Key features at a glance

- › Welded, all stainless steel construction
- › Edges hemmed for added rigidity and safety
- › 2" (5.1cm) ball bearing swivel casters
- › 16 gauge legs, 1.25" (3.18cm) outside diameter
- › Top and shelf stiffeners for additional strength
- › 3 year limited warranty

#### Specifications

Width: 72" (182.9cm)  
 Height: 34" (86.4cm)  
 Depth: 24" (61cm)  
 Distance between shelves: 19" (48.3cm)

#### Accessories

117 – Swivel casters locking 2" (5.1cm) (Set of 4)  
 480 – Rubber crutch tips in place of casters (table will be 2" (5.1cm)) lower