

Stainless Steel Instrument Table with Shelf

Model SS8018

UMF Medical stainless steel instrument tables are made with 20 gauge stainless steel shelves and 16 gauge stainless steel legs. All legs come with swivel casters and a 3-year limited warranty.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › Welded, all stainless steel construction
- › Edges hemmed for added rigidity and safety
- › 2" (5.1cm) ball bearing swivel casters
- › 16 gauge legs, 1.25" (3.18cm) outside diameter
- › Top and shelf stiffeners for additional strength
- › 3 year limited warranty

Specifications

Width: 20" (50.8cm)
 Height: 34" (86.4cm)
 Depth: 16" (40.6cm)
 Distance between shelves: 19" (48.3cm)

Accessories

117 – Swivel casters locking 2" (5.1cm) (Set of 4)
 165 – Stainless Steel push handle
 479 – Rubber crutch tips in place of casters (table will be 2" (5.1cm) lower)