

## Stainless Steel Mayo Stand with Foot Pedal

### Model SS8310

UMF Medical stainless steel Mayo Stand options feature a foot pedal height control and a tip-resistant weighted flat base that easily slides under O.R. tables. The units are mounted on swivel casters.

The SS8310 comes with a removeable 13" X 19" stainless steel tray.



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

### Key features at a glance

- › All stainless steel construction
- › Tip-resistant, weighted flat base
- › Removable stainless steel tray 13"x19" (33cm x 48.3cm)
- › 2" (5.1cm) ball bearing swivel casters
- › Foot pedal height control
- › 3 year limited warranty

### Specifications

Min. Height: 38.25" (97.2cm)  
Max. Height: 63.875" (162.2cm)  
Tray Dimensions: 13"x19" (33cm x 48.3cm)