

Stainless Steel Mayo Stand with Foot Pedal

Model SS8311

UMF Medical stainless steel Mayo Stand options feature a foot pedal height control and a tip-resistant weighted flat base that easily slides under O.R. tables. The units are mounted on swivel casters.

The SS8311 comes with a removeable 16" X 21" stainless steel tray.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › All stainless steel construction
- › Tip-resistant, weighted flat base
- › Removable stainless steel tray 16" x 21" (40.6cm x 53cm)
- › 2" (5.1cm) ball bearing swivel casters
- › Foot pedal height control
- › 3 year limited warranty

Specifications

Min. Height: 38.25" (97.2cm)
Max. Height: 63.875" (162.2cm)
Tray Dimensions: 16" x 21" (40.6cm x 53cm)