

## Stainless Steel I.V. Stand, Two Hook, Hand Wheel

### Model SS8343

The UMF Medical I.V. Pole comes with a telescopic stainless steel rod, secured in place by a positive-locking hand wheel. The I.V. poles have a 3-year limited warranty.

The SS8343 model has two hooks, with a 24" diameter, five-legged, non-tip plastic base equipped with swivel casters.



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

### Key features at a glance

- › Weight capacity: 15 lbs. each hook/35 lbs. total
- › 1" (2.53cm) diameter tubular stainless steel construction
- › 2 Hooks
- › 24" (61cm) non-tip aluminum base (5 legs)
- › 2" (5.1cm) swivel casters
- › Positive locking hand wheel
- › 3 year limited warranty

### Specifications

Extended height: 92.5" (235cm)  
Non-extended height: 50.5" (128.27cm)  
Base diameter: 24" (61cm)  
Caster diameter: 2" (5.1cm)  
Tube diameter: 1" (2.54cm)