



## Three Panel Privacy Screen

### Model 8424

UMF Medical's Three Panel Privacy Screen provides instant privacy with its sturdy frame and curtains. Four casters make it easy to transport and reposition. With three easy-to-clean panels and lightweight construction, this privacy screen can be carried and stored when not in use.



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

#### Key features at a glance

- › Folds for easy storage
- › Powder coat paint
- › Lightweight 3/4" steel tubing frame, 1" tubing legs for strength and stability
- › Curtain is constructed with a non-penetrable surface. Excellent tear strength; punctures will not produce running rips.
- › Minimal assembly required

#### Specifications

Width: 70"  
Folded width: 30"  
Unfolded width: 90"  
Weight: 40lbs

#### Bulk package includes

- › 10/pallet bundle
- › Completed frame
- › Completed legs with casters
- › Privacy curtains
- › Curtain rods
- › 2 bolts to attach legs to frame