

# Large Instrument/ Storage and Supply Cabinet

## Model 7142

The UMF Medical 7142 Large Instrument / Storage and Supply Cabinet has two lockable upper doors with glass panes, double drawers and double doors in the lower section, three adjustable glass shelves in the upper section, and one removable, non-adjustable metal shelf in the lower section. Comes with a 3-year limited warranty.



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

### Key features at a glance

- › 20 gauge steel construction
- › Removable, non-adjustable shelf in lower section
- › Heavy-duty concealed door hinges
- › Standard lock on upper right door
- › Safety release lever on upper left door
- › Adjustable leveling glides
- › 3 year limited warranty

### Specifications

Width: 31.5" (80cm)  
 Height: 66.5" (168.9cm)  
 Depth: 15.5" (39.4cm)  
 Net Weight: 172lb (78kg)  
 Shipping Weight: 248lb (112.5kg)  
 Shipping Dimensions: L48xW40xH77  
 (121.9cm x 101.6cm x 195.6cm)

### Accessories

070 – 110 Volt Light  
 071 – 220 Volt Light  
 502 – Double door lock (factory installed)  
 505 – Single drawer lock (factory installed)  
 507 – Additional keys (x2)