

## Two-Step Stainless Steel Step Stool

### Model SS8370

All UMF Medical foot stools boast high-quality, all-welded construction 16 gauge steel. Includes a non-slip corrugated rubber step tread, and no-tip offset stool legs with non-slip rubber tips. All foot stools come with a 3-year limited warranty and a non-corrosive option.



**ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA**

### Key features at a glance

- › High-quality, welded stainless steel construction
- › 16 gauge, 1" (2.54cm) outside diameter stainless steel tubing
- › 1.25" x .125" (3.2cm x .3cm) stainless steel band below step platform for rigidity
- › Non-slip rubber feet on stool legs
- › Non-slip rubber step treads
- › 3-year limited warranty

### Specifications

Width: 18" (45.7cm)  
Height: 15.63" (39.7cm)  
Depth: 25" (63.5cm)