

Stackable Stainless Steel Foot Stool

Model SS8380

All UMF Medical foot stools boast high-quality, all-welded construction 16 gauge steel. Includes a non-slip corrugated rubber step tread, and no-tip offset stool legs with non-slip rubber tips. All foot stools come with a 3-year limited warranty and a non-corrosive option.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › High-quality, welded stainless steel construction
- › Interlocking design
- › 16 gauge, 1" (2.54cm) outside diameter stainless steel tubing
- › 1.25" x .125" (3.2cm x .3cm) stainless steel band below step platform for rigidity
- › Non-slip rubber feet on stool legs
- › Non-slip rubber step treads
- › 3-year limited warranty

Specifications

Width: 18" (45.7cm)
 Height: 6" (15.2cm)
 Depth: 12" (30.5cm)

Shown with optional transport cart (UMF Medical item SS8381)

