

Two Function Bariatric Power Table

Model 5005


The UMF Medical 5005 Bariatric Treatment Table provides unmatched capacity and comfort in a stable, heavy duty all-steel design. The 5005 has an 800-lb. (362.9kg) patient weight capacity and a low 19" seat height. It also features a fully articulating power back section as well as power hi-lo control.



Colors



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA



River Rock



Sand Gray



Steel Blue



Midnight Blue

Key features at a glance

- › Ultra comfort seamless top with antimicrobial guard and built in pillow
- › ADA Compliant with 19" (48.3cm) low access
- › All-steel heavy duty construction
- › Power Hi-Lo
- › Two function hand control
- › Concealed paper roll holder 18" – 21" (45.7cm – 53.3cm)
- › Electro-mechanical actuator system
- › 10" (25.4cm) leg section
- › Large, full extension storage drawer on roller bearing drawer guides
- › 800lb (363kg) weight capacity
- › 3 year limited warranty

Specifications

Height: Min. 19" (48.3cm) at seat Max. 35" (88.9cm) at seat
 Width: 30" (76.2cm)
 Length: 70" (177.8cm) to 80" (203.2cm) with leg section extended
 Storage drawer: 16" D x 19" W x 4.5" H (40.6cm D x 48.3cm W x 11.4cm H)
 Electrical: 120 VAC, 60 HZ or 220/240 VAC, 60HZ
 Net Weight: 330lb (150kg)
 Shipping Weight: 441lb (200kg)
 Shipping Dimensions: L78xW36xH30 (198cm x 91.5cm x 76.2cm)

Accessories

422 – Foot Control, two function. Compatible w/ Hand Control, allowing both to be used simultaneously.
 802 – Paper cutter band
 810 – Restraint straps assembly
 856 – 220/240 Volt (factory installed)
 5005F – Fire retardant foam