

Stainless Steel Revolving Stool with Foot Ring

Model SS6704

UMF Medical's stools are made of high-quality, medical grade 304 stainless steel. All come standard with tubular footrest and non-slip material, flush mounted onto stool seat (can be manufactured without). All come with a limited 3-year warranty.



ALL UMF MEDICAL PRODUCTS ARE MADE IN THE USA

Key features at a glance

- › High-quality, welded stainless steel construction
- › Welded foot rest
- › Non-slip material mounted into recess of stool seat (flush finish)
- › Non-slip rubber feet on stool legs (device shown with optional casters)
- › Heavy duty threaded stem for seat height adjustment
- › 3-year limited warranty

Specifications

Seat diameter: 14" (35.6cm)
Foot rest dimensions: 16.25"x16.25" (41.3cm x 41.3cm)
Min. seat height: 19" (48.3cm)
Max. seat height: 25" (63.5cm)

Accessories

112 – Swivel Casters 2" (set of 4)
117 – Swivel Casters Locking 2" (set of 4)
SS6704NC – Non-corrosive steel

**Mild corrosion on the stem of the stool can take place if subject to moisture (SS6704 only)*